


Bisquert Open Model



$$Z = \sqrt{\zeta_m \chi_m} \coth \left(L \sqrt{\frac{\chi_m}{\zeta_m}} \right) \quad \text{where } \zeta_m = r_m \text{ and } \chi_m = r_k \parallel q_m$$

이 때 Constant Phase Element q_m 이

$$q_m = \frac{1}{Q_y s^{Q_a}} \quad \text{where } s = j\omega \text{ and } j = \sqrt{-1}$$

으로 표현되어지므로

$$\chi_m = \frac{r_k}{1 + r_k Q_y s^{Q_a}}$$

$$Z = \sqrt{\frac{r_m r_k}{1 + r_k Q_y s^{Q_a}}} \coth \left(L \sqrt{\frac{r_m}{r_k} (1 + r_k Q_y s^{Q_a})} \right)$$

이 된다.

이 모델을 Element A라 하고 ZMAN에 추가해보자. Appendix 1, 2, 3를 참고하면서 다음 테이블의 Steps을 따라 하자.

1. Add "A" element in SIMPLE category	
Category	SIMPLE
Model	A
Formula	$\coth(AL * \sqrt{Am / Ak * (1 + Ak * Ay * \text{pow}(s, Aa))}) * \sqrt{Am * Ak / (1 + Ak * Ay * \text{pow}(s, Aa))}$
Note	<p>Refer to Appendix 1 and 2</p> <ul style="list-style-type: none"> • $AL = L$ • $Am = r_m$ • $Ak = r_k$ • $Ay = Q_y$ • $Aa = Q_a$ • $s = j\omega$
Steps	<p>Refer to Fig. 1 ~ 3</p> <ol style="list-style-type: none"> (1) Open Model Editor dialog (2) Select SIMPLE item in the Category list (3) Select Model > New/Edit... in menu (4) Select A as Symbol (5) Type "Bisquert Open" in the Alias blank (6) Type above formula in the Formula blank (7) Click the OK button (8) Select Parameter tab and type default values
2. Make Rs-A-L1 Model in USER category	
Category	USER
Model	Rs-A-L1
Steps	<p>Refer to Fig 4 ~ 5</p> <ol style="list-style-type: none"> (1) Select USER item in the Category list (2) Select Model > New... in menu (3) Type "Rs-A-L1" in the Model blank (4) Click the OK button (5) Select Parameter tab and type default values




Fig 1




Fig 2




Fig 3




Fig 4




Fig 6

Appendix 1. List of Functions

Constants	
pi	Returns pi

Functions	
abs(x)	computes the absolute value of x
acos(x)	inverse cosine of x
acosh(x)	inverse hyperbolic cosine of x
acot(x)	inverse cotangent of x
acoth(x)	inverse hyperbolic cotangent of x
asin(x)	inverse sine of x
asinh(x)	inverse hyperbolic sine of x
atan(x)	inverse tangent of x
atan2(x,y)	inverse tangent of x/y
atanh(x)	inverse hyperbolic tangent of x
ceil(x)	computes the smallest integer greater than or equal to x
cos(x)	cosine of x
cosh(x)	hyperbolic cosine of x
cot(x)	cotangent of x
coth(x)	hyperbolic cotangent of x
csc(x)	cosecant of x
csch(x)	hyperbolic cosecant of x
deg(x)	converts radians to degrees ($x \cdot 180/\pi$)
e(x)	returns e or the argument multiplied by e
erf(x)	Error function
erfc(x)	Complementary Error function
exp(x)	e raised to the x power (exponential function)
factr(x)	factorial
floor(x)	computes the largest integer less than or equal to x
fract(x)	computes the fractional part of x
gamma(x)	Gamma function
gammai(a,x)	Incomplete Gamma function
getexp(x)	computes the exponent of a floating-point value
getman(x)	computes the mantissa of a floating-point value
int(x)	computes the integer part of x
ldexp(m,e)	computes a floating-point number from mantissa and exponent
ln(x)	natural logarithm of x (logarithm to the base e)

log(x,y)	logarithm of y to the base x
log10(x)	logarithm of x to the base 10
log2(x)	logarithm of x to the base 2
pi(x)	returns pi or the argument multiplied by pi
pow(x,y)	x raised to the y power
pow10(x)	10 raised to the x power
pow2(x)	2 raised to the x power
rad(x)	converts degrees to radians ($x \cdot \pi / 180$)
random(x,y)	generates random numbers within the specified range
sec(x)	secant of x
sech(x)	hyperbolic secant of x
sign(x)	returns the sign of x
sin(x)	sine of x
sinc(x)	$\sin(x)/x$
sinh(x)	hyperbolic sine of x
spike(x)	Spike function
sqrt(x)	computes the square root of x
square(x)	Square function
step(x)	Step function
tan(x)	tangent of x
tanh(x)	hyperbolic tangent of x

Appendix 2. Naming Convention for Circuit Element

Parameter naming convention for circuit element

ZMAN accepts only the following parameters:

A, AA, AB, AC, ... , Aa, Ab, Ac, ...
B, BA, BB, BC, ... , Ba, Bb, Bc, ...
...
Z, ZA, ZB, ZC, ..., Za, Zb, Zc, ...

- A parameter must consist of symbol and (or) an alphabetical character.
- A parameter name is case sensitive; i.e. AA and Aa are considered to be distinct parameters.

Appendix 3. Naming Convention in a model

Element naming convention in a model

Variables are distinguished by names. There are the following rules of naming variables:

- The first character should be a Symbol character(A to Z) and followed by alphabetical characters (capital and small letters), decimal digits, and underscore (_)
- A name is case sensitive; i.e. lowercase and uppercase letters are considered to be distinct characters
- There is no limit of name length; herewith all characters are significant