

# High Current Electrochemical Workstation ZIVE SHP1003



**±3V@100Amp High current application  
10uHz to 50kHz EIS capability**

*For  
Fuel Cells  
Supercapacitors  
Electrolysis*

# High Current Electrochemical Workstation SHP1003

The **ZIVE SHP1003** is a high current potentiostat/galvanostat/impedance analyzer with a maximum current of 100A. The **ZIVE SHP1003** is the best choice complete DC and impedance characterization of various energy source such as fuel cell, supercapacitor, electroplating, electrolysis, etc.

The system is designed under FPGA and DSP control with high speed capability.

## DAC Control

: Two sets of high speed 16bit DAC(50MHz) for offset & scanning & one set of 16 bit DAC(1MHz) for auxiliary analog output control

## ADC Reading

: Two sets of 16 bit 500kHz ADC for reading voltage/current and 4 channel 16 bit 250kHz ADCs for auxiliary data input such as temperature, auxiliary voltage etc. It provides high frequency EIS, fast pulse techniques and high speed sampling time.

The **ZIVE SHP1003** provides impedance measurements over the frequency range 10uHz to 50kHz. The system's major target application is for energy devices especially for Fuel cell application but it might be used in other electrochemical application requiring high current and low voltage . With various advanced software packages, user can widen **ZIVE SHP1003**'s flexibility.

## System Features

- Versatile high current Potentiostat/Galvanostat/Impedance Analyzer
- 16 EIS techniques capability including multisine & real time EIS at affordable price
- IR measurement in constant current control is available
- Main applications
  - Fuel cell, supercapacitor, electroplating, electrolysis, etc.
- User defined alias and unit display for auxiliary signal.
- High speed data sampling time
  - 2usec or 3usec depending on data point number
- 3 measurement/control voltage ranges &
- 9 measurement/control current ranges
- Internal 542,000 data point storage and continuing experiment regardless of PC failure.
- Full software packages are included as standard
  - EIS test software package(EIS)
  - Energy software package(BAT)
  - Corrosion test software package(COR)
  - Electrochemical analysis software package(EAS)
- Free software upgrade


## Hardware Features

- $\pm 3V@100A$  control range
- High current ranges(100A) for high current application
- Smart LCD display
- Independent operation by DSP with FPGA
- Simultaneous 3 auxiliary voltage measurements
- Temperature measurement as standard
- 1 auxiliary analog output
- 2 digital outputs & 1 digital input
- Separated power and sensing line
- Max 50kHz EIS capability

## • Front View



## • Channel View



## Versatility

The **ZIVE SHP1003** comes with additional 3 analog inputs (auxiliary voltage input) and 1 analog output along with 2 digital outputs and 1 digital input, and one temperature input for K type thermocouple. It will help users expand the usage of the instrument.

For example,

1. User can measure the voltage between working and reference electrode and, by using 3 additional analog inputs(auxiliary voltage input), user can also measure the voltage between reference and counter electrode and between working and counter electrode as well.
2. With analog output, the system can control rotating speed of a rotator, MFC flow rate etc. by  $\pm 3V$  full scale.
3. User can control on/off of max. 2 devices by DO etc.

## Safety and Maintenance

1. Even though the communication failure occurs between PC and ZIVE SHP1003 the system continues its experiment on channel and saves the data into ZIVE memory up to 542,000 data points set. After the communication is restored, ZIVE will transfer saved data to PC automatically or user can transfer data when he/she wants. This function will be highly efficient for long time experiment.
2. User can define a safety condition setting by inputting his/her own safety levels for voltage, current, temperature etc. If the measurement value exceeds this setting value, the system will automatically stop to protect the system and cell.
3. If the control value of voltage or current is different from measured value, the experiment will stop automatically to protect the cell.
4. Automatic calibration function is available for user calibration.
5. The system is controlled from a PC via USB.
6. There is an emergency button to cell off for emergency.

## Application

The ZIVE SHP1003 electrochemical workstation is the best choice to complete DC and impedance characterization of various energy source such as fuel cell, supercapacitor, electroplating, electrolysis, etc.

### Fuel Cells



The ZIVE SHP1003 is ideal for characterizing the fuel cells and anodic/cathodic process mechanism at development and research grade. This system can be directly used for PEMFC, DMFC, and DEFC etc. The FRA can control an external electronic load for EIS measurement of fuel cell. I-V curve measurements in a full range of available current (autorange option is active during the I-V scan in order to ensure measurement with continuously high resolution).

### Super Capacitors



The ZIVE SHP1003 has fast potentiostat circuit with high speed data acquisition (50 usec/point, burst mode). This function is well applicable to super capacitor testing. Charging/discharging capability is used for this application.

### Electrolysis



The ZIVE SHP1003 can be used for electrolysis experiments. Potentiostatic, Galvanostatic measurement can be used for this application.

### Electroplating



The ZIVE SHP1003 can be used for electroplating experiments. Potentiostatic, Galvanostatic, Pulse voltammetry, Pulse amperometry measurement can be used for this application.

## Main Software

The Smart Manager (SM) is to control ZIVE SHP1003 model and it provides user defined sequential test by using sequence file, technique menu and batch file. The batch file allows the users to do a serial test by combining sequence files and/or technique files.

The SM software is easy to use and supports various electrochemical experiments including functions of system control, schedule file editor, real time graph, analysis graph, user calibration, and data file treatment etc.




Technique list

## Basic Techniques

Basic techniques with standard functions

- 1) Potentiostatic
- 2) Galvanostatic
- 3) Double step potentiostatic
- 4) Double step galvanostatic
- 5) OCP measurement
- 6) Potential sweep
- 7) Current sweep
- 8) Cyclic voltammetry
- 9) Fast potential sweep
- 10) Potentiostatic Ru measurement
- 11) Galvanostatic Ru measurement



The above functions can be used sequentially by step control function.

# High Current Electrochemical Workstation SHP1003

## Sequence editor

User can design his/her experiment procedure by using TASK sequential routine editor.

## Control Task Parameters


Control Mode			
constant	GSTAT	constant current control	
	Crate	constant Crate control	
	PSTAT	constant voltage control	
	POWER	constant power control	
	LOAD	constant load control	
	CC-CV	constant current constant voltage control	
	Crate-CV	Crate constant voltage control	
	CP-CV	constant power constant voltage control	
	CL-CV	constant load constant voltage control	
	Id	Id control	
	Is	Is control	
	OCP	OCP control	
	Step	GSTAT	current step control
		PSTAT	potential step control
Sweep	GSTAT	current sweep control	
	FAST-G	fast current sweep control	
	PSTAT	potential sweep control	
	FAST-P	fast potential sweep control	
EIS	ACV	AC voltammetry	
	GSTAT	galvanostatic EIS	
	PSTAT	potentiostatic EIS	
	OCP	OCP EIS	
	PSEUDO	pseudo galvanostatic EIS	
	HFR G	galvanostatic HFR	
	HFR P	potentiostatic HFR	
	MsineG	galvanostatic multisine EIS	
MsineP	potentiostatic multisine EIS		
Pulse	RTIG	galvanostatic real time EIS	
	RTIP	potentiostatic real time EIS	
	Vpulse	voltage pulse control	
	Ipulse	current pulse control	
Rest	GSINE	current sine wave control	
	PSINE	potential sine wave control	
	ZRA	rest control	
Loop	ZRA	ZRA control	
	MUX	loop control	
Device	D OUT	Digital output control	
	A OUT	Analog output control	
	TEMP CTRL	Temperature controller control	
	MUX	Multiplexer control	

## Sampling Condition

- time, |dI/dt|, |dV/dt|, |dT/dt|, |dA1/dt|, burst time

## Flow View



- User can see the sequence flow at a glance.



Flow view

## Batch function

User can design batch file including multiple technique files and/or sequence files. With this batch file, user can experiment several techniques/sequence in series automatically.



Sequence editor

- Constant potential, current, C-rate, power, load, OCP
- Sweep potential, current
- Fast sweep potential, current
- Staircase potential, current
- CC-CV, CP-CV, CL-CV, Crate-CV control
- Id, Is control
- EIS control
- Pulse or sinewave control
- Rest(voltage monitoring only)
- Loop(cycle) control

## Cut-off(Vertex) Condition

- Time(step, test, loop, cycle)
- Current, current density
- Voltage
- Capacity
- C-rate
- dV
- |dV/dt|
- |dI/dt|
- Aux1
- Eoc
- etc.

Condition-1		
Item	OP	DeltaValue
Step Time	>=	30
None		
Step End		
Step Time		
Current		
I Density		
Voltage		
[Capacity]		
-dV		
dV/dt		
dI/dt		
dT/dt		
Temp.(C)		
Eoc		
WVH		
LCC(%)		
LCD(%)		
FCC(%)		
Power(W)		
SumQ(Ah)		
SumE(VWh)		
Loop Next		
DT Ch		
Crate		


Cutoff condition

## Smart Manager Advanced Software Package

For a wide range of application, advanced software packages for specific experimental techniques are available as standard. Each software package's upgrade will be provided at free of charge.


### EIS Software Package(EIS)

- Potentiostatic EIS
- Galvanostatic EIS
- Pseudo galvanostatic EIS
- OCP<sup>(\*)</sup> EIS
- Potentiodynamic PEIS
- Galvanodynamic GEIS
- Potentiodynamic HFR
- Galvanodynamic HFR
- Potentiostatic HFR monitor
- Galvanostatic HFR monitor
- Multisine potentiostatic EIS
- Multisine galvanostatic EIS
- Intermittent potentiostatic EIS
- Intermittent galvanostatic EIS
- Real time potentiostatic EIS
- Real time galvanostatic EIS



Potentiostatic EIS data of metal sample soaked in NaCl solution over time, 3D Nyquist plot by ZMAN


(\*) The system measures open circuit potential before for each frequency change and apply AC sine wave on this potential.




Galvanostatic EIS

Potentiodynamic PEIS

Rs, Cp & Idc vs Vdc plot



Potentiostatic EIS




Intermittent Potentiostatic EIS


## Corrosion Software Package(COR)

Corrosion technique supports IR compensation.


1. Tafel(Tafel experiment)
2. Rp(Polarization resistance)
3. Potentiodynamic
4. Galvanodynamic
5. Cyclic polarization
6. Ecorr vs. time
7. Galvanic corrosion
8. RpEc trend
9. Reactivation potential
10. Critical pitting potential
11. Critical pitting temperature(option)
12. ASTM critical pitting temperature(option)
13. Potentiostatic ECN
14. Galvanostatic ECN
15. ZRA mode ECN




Tafel experiment



Rp (Polarization resistance)




Cyclic polarization resistance




Galvanic Corrosion Test

## Electrochemical Analysis Software Package(EAS)


1. Step techniques
  - CA(Chronoamperometry)
  - CC(Chronocoulometry)
  - CP(Chronopotentiometry)
2. Sweep techniques
  - LSV(Linear sweep voltammetry)
  - SDV(Sampled DC voltammetry)
  - Fast CV
  - Fast LSV
3. Pulsed techniques
  - DPV(Differential pulse voltammetry)
  - SWV(Square wave voltammetry)
  - DPA(Diff. pulse amperometry)
  - NPV(Normal pulsed voltammetry)
  - RNPV(Reverse normal pulse voltammetry)
  - DNPV(Differential normal pulse voltammetry)



50usec sampling



Sampled DC voltammetry



NPV & RNPV overlay

## Battery Software Package(BAT)

BAT software supports IR measurement.

1. Battery test techniques
  - CC/CV test for cycle life test of lithium battery
  - CC/CC test for cycle life test of NiCd or NiMH battery
  - Discharging test
  - EVS(Electrochemical voltage spectroscopy)
  - Variable scan rate CV
  - Potentiostatic IV curve
  - Galvanostatic IV curve
  - Steady-state CV
  - GITT(Galvanostatic intermittent titration technique) test
  - PITT(Potentiostatic intermittent titration technique) test
  - Pulse mode is available for GSM & CDMA profile.  
Pulse shape profile can measured by user's demand.
2. Control mode
  - Charge: CC, CC-CV, pulse, sine wave
  - Discharge: CC, CP, CR, pulse, sine wave
3. Cutoff condition
  - time, voltage, current, power, temperature, auxV etc.

Various battery charge/discharge test is available including pulse discharge for GSM and CDMA application.


## Control & Real Time Graph

Smart Manager provides virtual control panel for control & data acquisition with real time graph.

User can control and monitor for specific channel in details and he/she can monitor data in VOI(value of interest) window and channel status in one window. Real time graph's X, Y axis format will be changed per technique automatically. It can be defined by user's demand per techniques.



For experiment using sequence file or batch file, user can designate X,Y parameter on three different real time graph. The real time graph's format can be also selected.

The real time graph and VOI will be changed depending on DC test or impedance test automatically. The virtual control panel always displays the graph for recent test result. For impedance measurement, wave monitor will be displayed on real time graph to check wave's quality. This monitor can be switched to Lissajous(I vs. E) plot.






### 1) DC Graph

- For general data display
- 4 shortcut buttons: I vs. V, E vs. LogI, V, I vs. time, V vs. Q
- Graph parameters: time, Eref, I, Eoc, Id, Aux1, Aux2, Aux3, temp, LogI, Load, ChQ, DchQ, ChQs, DchQs, Ch P, Dch P, Ch-Wh, Dch-Wh, Sum Wh, Sum Q, Sum |Q|, |Q|, Rp, dQ/dV




### 2) EIS Graph

- For EIS data display
- 3 shortcut buttons: Nyquist plot, Bode plot, Cs vs. frequency
- Graph parameters: Frequency, Zre, -Zim, Zmag, Zph, Y, Yimg, Y, |Y|, Yph, LogZ, LogY, Rs(R-C), Cs(R-C), Rp(R|C), Cp(R|C), Rs(R-L), Ls(R-L), Q(R-L), time, Vdc, Idc, temp, Aux(1,2,3)



3D Bode Plot by ZMAN  
Technique used: Potentiodynamic impedance measurement by using a corrosion cell

### Strip Chart





Strip chart recorder function provides real graph function independently. You can monitor 2 Y axis data such as voltage, current, auxV1,2,3, temperature, power, and capacity etc. in real time.

## Graph



Smart Manager's graph function is to simplify the operation. There are 3 kinds of graph per each experiment. You can change X, Y1, Y2, Y3, Y4 axis parameter as you want. Each graph provides shortcut buttons. When you click these buttons, the format of the graph will be changed accordingly.



In DC and Cycle graph, whenever you click or the parameters which are related to current such as current, capacity, energy, power, load, etc., are changed into calculated specific value or density value, respectively.


: value divided by weight

: value divided by active area

### 3) BAT Graph

- For battery cycle data display
- 3 shortcut buttons: cycle capacity, cycle average, Log(cycle No) vs. depth of discharge plot.
- Graph parameters: cycle number, Ch Q, Dch Q, Sum Q, Coulomb Eff, Ch-Wh, Dch-Wh, Sum Wh, Energy Eff, MinV, MaxV, ChQs, DchQ, ChVavg, DchVavg, Vavg

### Data Export to ASCII & Excel File




Selectable between 'Convert data on graph only' and 'Convert selected file(s)'

## Data Analysis Software


ZIVE data file can be used for analysis by using external IVMAN™ software for DC analysis, IVMAN DA™ software for battery data analysis, IVMAN PA™ software for photo-voltaic cell data analysis and ZMAN™ software for EIS data analysis without license.

### ZMAN™ EIS Data Analysis Software


- Model simulation and fitting
- 2D- and 3D-Bode- and Nyquist plots
- Automatic equivalent circuit model search function
- Project concept to handle multiple EIS data analysis
- Parameter plot from fitted elements value
- Compatible with data format from Zahner, Gamry, Ametek etc. (License code is needed.)
- Various weighting algorithm
- Model library and user model
- KK plot
- Batch fitting for project data
- Impedance parameter simulation
- Interpolate bad data
- Black-Nichols plot
- 3D graph setting option
- Improved model editor
- Application model library for automatic searching
- Parameter simulation of model
- Genetic algorithm option for initial guessing
- Automatic initial guessing
- Trace movie function on fitting
- Free for ZIVE's data format(\*.seo, \*.wis) analysis (No license code required.)
- Circle fitting
- Data editing available (insert, delete, edit)
- Add/subtract element parameters
- Add/subtract model parameters
- Impedance, Z in polar, admittance, Y in Polar, modulus, M in polar, dielectric constant, E in polar. data display
- Empty cell capacitance calculation
- Find file function
- Data replacement by formula function
- Cursor data display
- Model finding result automatic sorting by Chi square value
- R, C R, L R, Q preview & graphic
- ZHIT function
- Mott-Schottky analysis
- Donor density vs. Vfb graph
- C vs. voltage graph





Cursor data display




Data replacement by formula function




Fitting display




Model editor & model library




Automatic model searching





LEVM fitting




Parameter simulation




Element add/subtraction




Parameter plot




Finding data file menu




Circular fitting




Importing 3rd parties ASCII data file



Project manager with data preview




2D Nyquist plot




2D Bode plot

# High Current Electrochemical Workstation SHP1003





Rp, Cp vs frequency (R | C)




Empty cell capacitance

• Tafel calculation result





Donor density vs. Vfb graph and analysis




KK consistency


• Tafel region selection & data list



Mott-Schottky analysis window



Bode & Nyquist overlay & 3D plots



C/R-V graph




## IVMAN DA™ Battery Test Data Analysis Software

- Battery test data analysis
- Electrochemical voltage spectroscopy (dQ/dV vs. V)
- Voltage vs. Capacity analysis (V vs. Q)
- Cycle graph (Q vs. cycle)
- Differential voltage graph (dV/dQ vs. Q)



## IVMAN™ Photovoltaic Cell Analysis



- Automatic analysis of parameters
- open circuit voltage, open circuit current, max. power, efficiency
- photo induced current, diode quality factor, series resistance, etc.



## IVMAN™ Main Software

- Ideal for DC corrosion data analysis and electro-analytical data analysis
- Initial guessing function on Tafel analysis
- Polarization resistance fitting
- 3D graph
- Find peak function
- Interpolation, differentiation, integration etc.
- Reporting function

## IVMAN™ DC Data Analysis Software





IVMAN™ software package consists of

- IVMAN software
- IVMAN utilities
- IVMAN main software
- IVMAN differential analysis software
- IVMAN photovoltaic cell analysis
- IVMAN Tafel analysis
- IVMAN extractor
- IVMAN peak find module




## IVMAN TA™ Tafel Analysis

- Simple Tafel calculation




Polarization resistance fitting



Polarization analysis result


## Optional Accessories

### • Single Cell Hardware Fixture




- for PEMFC & DMFC
- max. temp. : 120 °C or 180 °C
- active area : 5, 25cm<sup>2</sup>
- MEA is not included.


### • Membrane Conductivity Cell




- for 5, 25cm<sup>2</sup> fuel cell hardware fixture
- material : PEEK(cell body), platinum(wire)
- operating temp. : up to 130 °C




Tafel analysis




Time graph



Find peak menu




CV graph



3D graph




Edit data menu



Universal graph


## IVMAN EX™ Extractor

- Extracting data by cycle number or step
- Exporting ASCII file



## IVMAN PF™ Peak Find Module

- Independent peak finding software



## • Jig

- 4 pin probe knob type
- lever or knob type

**4PL-CBJ8-8P**



- 4 pin lever type cylindrical cell jig
- 8 channels, Max.current 50A
- Direct connection method for jig cables
- For high current application

**UCBJ1**



- 4 pin lever type high current universal cell jig
- For 1 channel, Max.current 50A
- Banana connectors for voltage measurement and connectors for high current

**4PL-UCBJ4-4P**



- 4 pin lever type high current universal cell jig
- For 4 channels, Max.current 50A
- 4 channels per panel
- For high current application

**4PK-UCJH4-4P**



- 4 pin knob type high current universal cell jig
- For 4 channels, Max.current 50A
- 4 channels per panel
- For high current application



High current cylindrical battery holder, P/N: HCCBJ65L



High current cylindrical battery holder, P/N: HCCBJ100L

## Specification

Main System	
PC communication	USB2.0 high speed
Line voltage	100~240VAC, 50/60Hz
Power requirement	1000Watt
Size	447.1x600x241mm(WxDxH), 28kg
LED indicator	Run, Comm

System	
Cell cable set	1.5 meter power cables (WE/CE CAMLOK terminal) & Reference, working sense, AuxV1 cable
Control DAC	DSP with FPGA 2x16bit DAC(50MHz) for bias & scan 1x16bit DAC(1MHz) for analog output
Data acquisition ADC	2x16bit ADCs(500kHz) for voltage, current 4x16bit ADCs(250kHz) for auxiliary voltage and temperature reading
Calibration	Automatic
Filter selection	4ea (5Hz, 1kHz, 500kHz, 5MHz)
Scan rate	0~200V/sec in common mode 0~5000V/sec in fast mode
Internal data memory	542,000 points
LCD display	DC & EIS mode automatically

Power Amplifier(CE)	
Compliance voltage	±3V
Max. current	±100A
Control speed selection	4ea
Bandwidth	200kHz
Slew rate	2V/μsec

Potentiostat Mode (voltage control)	
Voltage control	
Control voltage range	±3V, ±300mV, ±30mV
Voltage resolution	16 bit per each range
Voltage accuracy	±0.03% f.s. (gain x1)
Max. scan range	±3V vs. ref. E
Current measurement	
Current range	10 ranges 100A~100nA 100nA with gain
Current resolution	16 bit 3mA, 300μA, 30μA, 3μA, 300nA, 30nA, 3nA, 300pA, 30pA, 3pA
Current accuracy	±0.05% f.s.(gain x1)> 10uAf.s

Galvanostat Mode (current control)	
Current control	
Control current range	max. ±100A ± full scale depending on selected range
Current resolution	16 bit 3mA, 300μA, 30μA, 3μA, 300nA, 30nA, 3nA, 300pA, 30pA, 3pA
Current accuracy	±0.05% f.s.(gain x1)> 10uA f.s.
Voltage measurement	
Voltage range	±3V, ±300mV, ±30mV
Voltage resolution	16 bit 90uV, 9uV, 900nV
Voltage accuracy	±0.03% f.s. (gain x1)

Electrometer	
Max. input voltage	±10V
Input impedance	>2x10 <sup>13</sup> Ω   4.5pF
Bandwidth	>22MHz
CMRR	>114dB

EIS Measurement for System	
Frequency range	10uHz~50kHz
Frequency accuracy	<0.01%
Frequency interval setting (Point/Decade)	Max. 1000/decade (<43mHz) Max. 5000/decade (>43mHz)
Amplitude	0.5mV ~ 2.1V rms (Potentiostatic) 0.1~70% f.s. (Galvanostatic)
Mode	<b>Static EIS:</b> Potentiostatic, Galvanostatic, Pseudogalvanostatic, OCP <b>Dynamic EIS:</b> Potentiodynamic, Galvanodynamic <b>Fixed frequency impedance:</b> Potentiostatic, Galvanostatic, Potentiodynamic, Galvanodynamic <b>Multisine EIS:</b> Potentiostatic, Galvanostatic <b>Intermittent PEIS/GEIS</b>

Interfaces for System	
Digital & Auxiliary port	
Digital output	3 (open collector)
Digital input	1 (photo coupler)
Auxiliary voltage inputs	3 analog inputs: ±10V For measurement of WE vs. CE CE vs. RE or other signal
Analog output	1 analog output: ±10V For stirrer, MFC, RDE, etc.
Peripheral communication	I2C to control external devices
Temperature input port	
Temp. measurement	1 K-type thermocouple input

Software	
Max. step per experiment	1000
Shutdown safety limits	Voltage, Current, Temperature, etc.
Max. sampling rate	2μsec or 3μsec depending on data point number
Min. sampling time	Unlimited
Sampling condition	Time, dv/dt, dl/dt, temperature, etc.

PC Requirement	
Operating system	Windows 7/8/10(32bit/64bit OS)
PC specification	Pentium4, RAM 1GB or higher
Display	1600x900 high color or higher
USB	High speed 2.0

General	
Dummy cell	One external dummy cell included
Thermocouple	K-type, 1.5 meter long (option)
Auxiliary cable	Option
Misc. cable	Option
Impedance analysis S/W	ZMAN™ software
DC data analysis S/W	IVMAN™ software package

The specifications are subject to change without notice.  
Windows is a registered trademark of Microsoft Corporation.

Designed by

**ZIVE LAB**  
www.zivelab.com

**WonATech**

WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-2-576-2635  
e-mail) sales@wonatech.com  
website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified