

High Power Electrochemical Workstation ZIVE SHP1003



$\pm 3\text{V}@100\text{Amp}$ power application
10uHz to 50kHz EIS capability

*For
Fuel Cell
Supercapacitor
Electroplating
Electrolysis*

Electrochemical Workstation ZIVE SHP1003

The **ZIVE SHP1003** is a high current potentiostat/galvanostat/impedance analyzer with a maximum current of 100A. The **ZIVE SHP1003** is the best choice complete DC and impedance characterization of various energy source such as fuel cell, supercapacitor, electroplating, electrolysis, etc.

The system is designed under FPGA and DSP control with high speed capability.

DAC Control

: Two sets of high speed 16bit DAC(50MHz) for offset & scanning & one set of 16 bit DAC(1MHz) for auxiliary analog output control

ADC Reading

: Two sets of 16 bit 500kHz ADC for reading voltage/current and 4 channel 16 bit 250kHz ADCs for auxiliary data input such as temperature, auxiliary voltage etc. It provides high frequency EIS, fast pulse techniques and high speed sampling time.

The **ZIVE SHP1003** provides impedance measurements over the frequency range 10uHz to 50kHz. The system's major target application is for energy devices especially for Fuel cell application but it might be used in other electrochemical application requiring high current low voltage. With various advanced software packages, user can widen **ZIVE SHP1003**'s flexibility.

System Features

- Versatile high current Potentiostat/Galvanostat/Impedance Analyzer
- 18 EIS techniques capability including multisine & real time EIS at affordable price
- IR measurement in constant current control is available
- Main applications
 - Fuel cell, supercapacitor, electroplating, electrolysis, etc.
- User defined alias and unit display for auxiliary signal.
- High speed data sampling time
 - 2usec or 3usec depending on data point number
- 3 measurement/control voltage ranges &
- 9 measurement/control current ranges
- Internal 542,000 data point storage and continuing experiment regardless of PC failure.
- Full software packages are included as standard
 - EIS test software package(EIS)
 - Energy software package(BAT)
 - Corrosion test software package(COR)
 - Electrochemical analysis software package(EAS)
- Free software upgrade


Hardware Features

- $\pm 3V@100A$ control range
- High current ranges(100A) for high current application
- Smart LCD display
- Independent operation by DSP with FPGA
- Simultaneous 3 auxiliary voltage measurements
- Temperature measurement as standard
- 1 auxiliary analog output
- 2 digital outputs & 1 digital inputs
- Separated power and sensing line
- Max 50kHz EIS capability

• Front View



• Channel View



Versatility

The **ZIVE SHP1003** comes with additional 3 analog inputs (auxiliary voltage input) and 1 analog output along with 2 digital outputs and 1 digital input, and one temperature input for K type thermocouple. It will help users expand the usage of the instrument.

For example,

1. User can measure the voltage between working and reference electrode and, by using 3 additional analog inputs(auxiliary voltage input), user can also measure the voltage between reference and counter electrode and between working and counter electrode as well.
2. With analog output, the system can control rotating speed of a rotator, MFC flow rate etc. by $\pm 3V$ full scale.
3. User can control on/off of max. 2 devices by DO etc.
4. This system can interface with an internal booster.

Safety and Maintenance

1. Even though the communication failure occurs between PC and ZIVE SHP1003 the system continues its experiment on channel and saves the data into ZIVE memory up to 542,000 data point set. After the communication is restored, ZIVE will transfer saved data to PC automatically or user can transfer data when he/she wants. This function will be highly efficient for long time experiment.
2. User can define a safety condition setting by inputting his/her own safety levels for voltage, current, temperature etc. If the measurement value exceeds this setting value, the system will automatically stop to protect the system and cell.
3. If the control value of voltage or current is different from measured value, the experiment will stop automatically to protect the cell.
4. Automatic calibration function is available for user calibration.
5. The system is controlled from a PC via USB.
6. There is an emergency button to cell off for emergency.

Application

The ZIVE SHP1003 electrochemical workstation is the best choice to best choice complete DC and impedance characterization of various energy source such as fuel cell, supercapacitor, electroplating, electrolysis, etc.

Fuel Cells



The ZIVE SHP1003 is ideal for characterizing the fuel cells and anodic/cathodic process mechanism at development and research grade. This system can be directly used for PEMFC, DMFC, and DEFC etc. The FRA can control an external electronic load for EIS measurement of fuel cell. I-V curve measurements in a full range of available current (autorange option is active during the I-V scan in order to ensure measurement with continuously high resolution).

Super Capacitors



The ZIVE SHP1003 has fast potentiostat circuit with high speed data acquisition (50usec/point, burst mode). This function is well applicable to super capacitor testing. Charging/discharging capability is used for this application.

Electroplating



The ZIVE SHP1003 can be used for electroplating experiments.

Potentiostatic, Galvanostatic, Pulse voltammetry, Pulse Amperometry measurement can be used for this application.

Electrolysis




The ZIVE SHP1003 can be used for electrolysis experiments.

Potentiostatic, Galvanostatic measurement can be used for this application.

Main Software

The Smart Manager (SM) is to control ZIVE SHP1003 model and it provides user defined sequential test by using sequence file, technique menu and batch file. The batch file allows the users to do a serial test by combining sequence files and/or technique files.

The SM software is easy to use and supports various electrochemical experiments including functions of system control, schedule file editor, real time graph, analysis graph, user calibration, and data file treatment etc.




Technique list

Basic Techniques

Basic techniques with standard functions

- 1) Potentiostatic
- 2) Galvanostatic
- 3) Double step potentiostatic
- 4) Double step galvanostatic
- 5) OCP measurement
- 6) Potential sweep
- 7) Current sweep
- 8) Cyclic voltammetry
- 9) Fast potential sweep
- 10) Potentiostatic Ru measurement
- 11) Galvanostatic Ru measurement




The above functions can be used sequentially by step control function.

Sequence editor

User can design his/her experiment procedure by using TASK sequential routine editor.

Control Task Parameters

| Control Mode | | |
|--------------|-----------|---|
| constant | GSTAT | constant current control |
| | Crate | constant Crate control |
| | PSTAT | constant voltage control |
| | POWER | constant power control |
| | LOAD | constant load control |
| | CC-CV | constant current constant voltage control |
| | Crate-CV | Crate constant voltage control |
| | CP-CV | constant power constant voltage control |
| | CL-CV | constant load constant voltage control |
| | Id | Id control |
| Step | Is | Is control |
| | OCP | OCP control |
| | GSTAT | current step control |
| Sweep | PSTAT | potential step control |
| | GSTAT | current sweep control |
| | FAST-G | fast current sweep control |
| EIS | PSTAT | potential sweep control |
| | FAST-P | fast potential sweep control |
| | ACV | AC voltammetry |
| Pulse | GSTAT | galvanostatic EIS |
| | PSTAT | potentiostatic EIS |
| | OCP | OCP EIS |
| | PSUEDO | pseudo galvanostatic EIS |
| | HFR G | galvanostatic HFR |
| | HFR P | potentiostatic HFR |
| | M sineG | galvanostatic multisine EIS |
| | M sineP | potentiostatic multisine EIS |
| | RTIG | galvanostatic real time EIS |
| | RTIP | potentiostatic real time EIS |
| Rest | Vpulse | voltage pulse control |
| | Ipulse | current pulse control |
| | GSINE | current sine wave control |
| | PSINE | potential sine wave control |
| ZRA | | rest control |
| | ZRA | ZRA control |
| Loop | | loop control |
| | D OUT | Digital output control |
| | A OUT | Analog output control |
| Device | TEMP CTRL | Temperature controller control |
| | MUX | Multiplexer control |




Sequence editor

Sampling Condition

- time, |dI/dt|, |dV/dt|, |dT/dt|, |dA1/dt|, burst time

Flow View


- User can see the sequence flow at a glance.



Flow view

Batch function

User can design batch file including multiple technique files and/or sequence files. With this batch file, user can experiment several techniques/sequence in series automatically.




Smart Manager Advanced Software Package

For a wide range of application, advanced software packages for specific experimental techniques are available as standard. Each software package's upgrade will be provided at free of charge.

EIS Software Package(EIS)

- Potentiostatic EIS
- Galvanostatic EIS
- Pseudo galvanostatic EIS
- OCP(*) EIS
- Potentiodynamic PEIS
- Galvanodynamic GEIS
- Potentiodynamic HFR
- Galvanodynamic HFR
- Potentiostatic HFR monitor
- Galvanostatic HFR monitor
- Multisine potentiostatic EIS
- Multisine galvanostatic EIS
- Intermittent potentiostatic EIS
- Intermittent galvanostatic EIS
- Real time potentiostatic EIS
- Real time galvanostatic EIS

(*) The system measures open circuit potential before for each frequency change and apply AC sine wave on this potential.



Potentiostatic EIS data of metal sample soaked in NaCl solution over time, 3D Nyquist plot by ZMAN


- Constant potential, current, C-rate, power, load, OCP
- Sweep potential, current
- Fast sweep potential, current
- Staircase potential, current
- CC-CV, CP-CV, CL-CV, Crate-CV control
- Id, Is control
- EIS control
- Pulse or sinewave control
- Rest(voltage monitoring only)
- Loop(cycle) control

Cut-off(Vertex) Condition

- Time(step, test, loop, cycle)
- Current, current density
- Voltage
- Capacity
- C-rate
- dV
- |dV/dt|
- |dI/dt|
- Aux1
- Eoc
- etc.

| Condition-1 | | |
|-------------|----|------------|
| Item | OP | DeltaValue |
| Step Time | >= | 30 |
| None | | |
| Step End | | |
| Step Time | | |
| Current | | |
| I Density | | |
| Voltage | | |
| Capacity | | |
| -dV | | |
| dV/dt | | |
| dI/dt | | |
| Temp.(C) | | |
| AUX1 | | |
| AUX2 | | |
| AUX3 | | |
| Test Time | | |
| Loop Time | | |
| Cycle Time | | |
| Eoc | | |
| WHr | | |
| LCC(%) | | |
| FCC(%) | | |
| PCC(%) | | |
| Power(W) | | |
| SumQ(Ah) | | |
| SumE(WWh) | | |
| Loop Next | | |
| OT Ch | | |
| Crate | | |


Cutoff condition




Galvanostatic EIS

Potentiodynamic PEIS

Rs, Cp & Idc vs Vdc plot



Potentiostatic EIS




Intermittent Potentiostatic EIS


■ Corrosion Software Package(COR)

Corrosion technique supports IR compensation.


1. Tafel(Tafel experiment)
2. Rp(Polarization resistance)
3. Potentiodynamic
4. Galvanodynamic
5. Cyclic polarization
6. Ecorr vs. time
7. Galvanic corrosion
8. RpEc trend
9. Reactivation potential
10. Critical pitting potential
11. Critical pitting temperature(option)
12. ASTM critical pitting temperature(option)
13. Potentiostatic ECN
14. Galvanostatic ECN
15. ZRA mode ECN




Tafel experiment



Rp (Polarization resistance)




Cyclic polarization resistance




Galvanic Corrosion Test

■ Electrochemical Analysis Software Package(EAS)


1. Step techniques
 - CA(Chronoamperometry)
 - CC(Chronocoulometry)
 - CP(Chronopotentiometry)
2. Sweep techniques
 - LSV(Linear sweep voltammetry)
 - SDV(Sampled DC voltammetry)
 - Fast CV
 - Fast LSV
3. Pulsed techniques
 - DPV(Differential pulse voltammetry)
 - SWV(Square wave voltammetry)
 - DPA(Diff. pulse amperometry)
 - NPV(Normal pulsed voltammetry)
 - RNPV(Reverse normal pulse voltammetry)
 - DNPV(Differential normal pulse voltammetry)



50usec sampling



Sampled DC voltammetry



NPV & RNPV overlay

■ Battery Software Package(BAT)

BAT software supports IR measurement.

1. Battery test techniques
 - CC/CV test for cycle life test of lithium battery
 - CC/CC test for cycle life test of NiCd or NiMH battery
 - Discharging test
 - EVS(Electrochemical voltage spectroscopy)
 - Variable scan rate CV
 - Potentiostatic IV curve
 - Galvanostatic IV curve
 - Steadystate CV
 - GITT(Galvanostatic intermittent titration technique) test
 - PITT(Potentiostatic intermittent titration technique) test
 - Pulse mode is available for GSM & CDMA profile.
Pulse shape profile can be measured by user's demand.
2. Control mode
 - Charge: CC, CC-CV, pulse, sine wave
 - Discharge: CC, CP, CR, pulse, sine wave
3. Cutoff condition
 - time, voltage, current, power, temperature, auxV etc.

Various battery charge/discharge test is available including pulse discharge for GSM and CDMA application.


Control & Real Time Graph

Smart Manager provides virtual control panel for control & data acquisition with real time graph.


User can control and monitor for specific channel in details and he/she can monitor data in VOI(value of interest) window and channel status in one window. Real time graph's X,Y axis format will be changed per technique automatically. It can be defined by user's demand per techniques.

For experiment using sequence file or batch file, user can designate X,Y parameter on three different real time graph. The real time graph's format can be also selected.

The real time graph and VOI will be changed depending on DC test or impedance test automatically. The virtual control panel always displays the graph for recent test result. For impedance measurement, wave monitor will be displayed on real time graph to check wave's quality. This monitor can be switched to Lissajous(I vs. E) plot.



Strip Chart





Strip chart recorder function provides real graph function independently. You can monitor 2 Y axis data such as voltage, current, auxV1,2,3, temperature, power, and capacity etc. in real time.

Graph



Smart Manager's graph function is to simplify the operation. There are 3 kinds of graph per each experiment. You can change X, Y1, Y2, Y3, Y4 axis parameter as you want. Each graph provides shortcut buttons. When you click these buttons, the format of the graph will be changed accordingly.





In DC and Cycle graph, whenever you click or the parameters which are related to current such as current, capacity, energy, power, load, etc., are changed into calculated specific value or density value, respectively.

: value divided by weight

: value divided by active area




1) DC Graph

- For general data display
- 4 shortcut buttons: I vs. V, E vs. LogI, V, I vs. time, V vs. Q
- Graph parameters: time, Eref, I, Eoc, Id, Aux1, Aux2, Aux3, temp, LogI, Load, ChQ, DchQ, ChQs, DchQs, Ch P, Dch P, Ch-Wh, Dch-Wh, Sum Wh, Sum Q, Sum |Q|, |Q|, Rp, dQ/dV



2) EIS Graph

- For EIS data display
- 3 shortcut buttons: Nyquist plot, Bode plot, Cs vs. frequency
- Graph parameters: Frequency, Zre, -Zim, Zmag, Zph, Y, Yimg, Y, |Y|, Yph, LogZ, LogY, Rs(R-C), Cs(R-C), Rp(R|C), Cp(R|C), Rs(R-L), Ls(R-L), Q(R-L), time, Vdc, Idc, temp, Aux(1,2,3)




3D Bode Plot by ZMAN
Technique used: Potentiodynamic impedance measurement by using a corrosion cell

3) BAT Graph

- For battery cycle data display
- 3 shortcut buttons: cycle capacity, cycle average, Log(cycle No) vs. depth of discharge plot.
- Graph parameters: cycle number, Ch Q, Dch Q, Sum Q, Coulomb Eff, Ch-Wh, Dch-Wh, Sum Wh, Energy Eff, MinV, MaxV, ChQs, DchQ, ChVavg, DchVavg, Vavg

Data Export to ASCII & Excel File





Selectable between 'Convert data on graph only' and 'Convert selected file(s)'

Data Analysis Software


ZIVE data file can be used for analysis by using external IVMAN™ software for DC analysis, IVMAN DA™ software for battery data analysis, IVMAN PA™ software for photo-voltaic cell data analysis and ZMAN™ software for EIS data analysis without license.

ZMAN™ EIS Data Analysis Software


- Model simulation and fitting
- 2D- and 3D-Bode- and Nyquist plots
- Automatic equivalent circuit model search function
- Project concept to handle multiple EIS data analysis
- Parameter plot from fitted elements value
- Compatible with data format from Zahner, Gamry, Ametek etc. (License code is needed.)
- Various weighting algorithm
- Model library and user model
- KK plot
- Batch fitting for project data
- Impedance parameter simulation
- Interpolate bad data
- Black-Nichols plot
- 3D graph setting option
- Improved model editor
- Application model library for automatic searching
- Parameter simulation of model
- Genetic algorithm option for initial guessing
- Automatic initial guessing
- Trace movie function on fitting
- Free for ZIVE's data format(*.seo, *.wis) analysis (No license code required.)
- Circle fitting
- Data editing available (insert, delete, edit)
- Add/subtract element parameters
- Add/subtract model parameters
- Impedance, Z in polar, admittance, Y in Polar, modulus, M in polar, dielectric constant, E in polar. data display
- Empty cell capacitance calculation
- Find file function
- Data replacement by formula function
- Cursor data display
- Model finding result automatic sorting by Chi square value
- R, C, L, R, Q preview & graphic
- ZHIT function
- Mott-Schottky analysis
- Donor density vs. Vfb graph
- C vs. voltage graph




Parameter plot




Importing 3rd parties ASCII data file




Project manager with data preview





Cursor data display




Data replacement by formula function




Fitting display




Model editor & model library




Automatic model searching




LEVM fitting




Parameter simulation




Element add/subtraction




Finding data file menu



Circular fitting




2D Nyquist plot




2D Bode plot

Electrochemical Workstation ZIVE SHP1003





Rp,Cp vs frequency (R | C)




Empty cell capacitance

- Tafel calculation Result





Donor density vs. Vfb graph and analysis




KK consistency


- Tafel region selection & data list



Mott-Schottky analysis window



Bode & Nyquist overlay & 3D plots



C/R-V graph

IVMAN™ DC Data Analysis Software




IVMAN™ software package consists of

- IVMAN software
- IVMAN utilities
 - IVMAN main software
 - IVMAN differential analysis software
 - IVMAN photo voltaic cell analysis.
 - IVMAN Tafel analysis
 - IVMAN extractor
 - IVMAN peak find module



IVMAN TA™ Tafel Analysis

- Simple Tafel calculation




IVMAN DA™ Battery Test Data Analysis Software

- Battery test data analysis
- Electrochemical voltage spectroscopy (dQ/dV vs. V)
- Voltage vs. Capacity analysis (V vs. Q)
- Cycle graph (Q vs. cycle)
- Differential voltage graph (dV/dQ vs. Q)



IVMAN™ Photovoltaic Cell Analysis




- Automatic analysis of parameters
 - open circuit voltage, open circuit current, max. power, efficiency
 - photo induced current, diode quality factor, series resistance, etc.




IVMAN™ Main Software

- Ideal for DC corrosion data analysis and electro-analytical data analysis
- Initial guessing function on Tafel analysis
- Polarization resistance fitting
- 3D graph
- Find peak function
- Interpolation, differentiation, integration etc.
- Reporting function




Polarization resistance fitting



Polarization analysis result


Optional Accessories

• Single Cell Hardware Fixture



- for PEMFC & DMFC
- max. temp. : 120 °C or 180 °C
- active area : 5, 25cm²
- MEA is not included.

• Membrane Conductivity Cell




- for 5, 25cm² fuel cell hardware fixture
- material : PEEK(cell body), platinum(wire)
- operating temp. : up to 130 °C


• Through-Plane Conductivity Test Jig



- for through plane conductivity measurement
- 2 probe type
- Sample size : >30mm diameter
- Sample thickness : max. 40mm




Tafel analysis




Time graph

Find peak menu



CV graph

3D graph




Edit data menu

Universal graph


IVMAN EX™ Extractor

- Extracting data by cycle number or step
- Exporting ASCII file



IVMAN PF™ Peak Find Module


- Independent peak finding software




Electrochemical Workstation ZIVE SHP1003

• Battery&Coin Cell Jig

- 4 pin probe knob type
- lever or knob type



• Cell Kit



• Universal Electrode Holder



- electrode and glass vial are not included.

• Faraday Cage



- size : 300 x 300 x 398mm (WxDxH)

Specification

| Main System | |
|-------------------|------------------------------|
| PC communication | USB2.0 high speed |
| Line voltage | 100~240VAC, 50/60Hz |
| Max. output power | 300Watt |
| Size | 447.1x600x241mm(WxDxH), 28kg |
| LED indicator | Run, Comm |

| System | |
|----------------------|---|
| Cell cable | 1 meter shielded type(standard) working, reference, counter, working sense |
| Control DAC | DSP with FPGA 2x16bit DAC(50MHz) for bias & scan 1X16bit DAC(1MHz) for analog output |
| Data acquisition ADC | 2x16bit ADCs(500kHz) for voltage, current 4x16bit ADCs(250kHz) for auxiliary voltage and temperature reading |
| Calibration | Automatic |
| Filter selection | 4ea(5Hz, 1kHz, 500kHz, 5MHz) |
| Scan rate | 0~200V/sec in common mode 0~5000V/sec in fast mode |
| Internal data memory | 542,000 points |
| LCD display | DC & EIS mode automatically |

| Power Amplifier(CE) | |
|-------------------------|------------------|
| Power | 300Watt(3V@100A) |
| Compliance voltage | ±3V |
| Max. current | ±100A |
| Control speed selection | 4ea |
| Bandwidth | 1MHz |
| Slew rate | V/usec |

| Potentiostat Mode (voltage control) | |
|-------------------------------------|--|
| Voltage control | |
| Control voltage range | ±3V, ±300mV, ±30mV |
| Voltage resolution | 16 bit per each range |
| Voltage accuracy | ±0.02% fs (gain x1) |
| Max. scan range | ±3V vs. ref. E |
| Current measurement | |
| Current range | 10 ranges 100A~100nA 100nA with gain |
| Current resolution | 16 bit 3mA, 300uA, 30uA, 3uA, 300nA, 30nA, 3nA, 300pA, 30pA, 3pA |
| Current accuracy | ±0.05%fs(gainx1)>10uAf.s |

| Galvanostat Mode (current control) | |
|------------------------------------|--|
| Current control | |
| Control current range | max. ±100A ± full scale depending on selected range |
| Current resolution | 16 bit 3mA, 300uA, 30uA, 3uA, 300nA, 30nA, 3nA, 300pA, 30pA, 3pA |
| Current accuracy | ±0.05%fs(gainx1)>10uAf.s |
| Voltage measurement | |
| Voltage range | ±3V, ±300mV, ±30mV |
| Voltage resolution | 16 bit 90uV, 9uV, 900nV |
| Voltage accuracy | ±0.02% fs(gain x1) |

| Electrometer | |
|--------------------|-------------------------------|
| Max. input voltage | ±10V |
| Input impedance | >2x10 ¹³ Ω 4.5pF |
| Bandwidth | >22MHz |
| CMRR | >114dB |

| EIS Measurement for System | |
|----------------------------|---|
| Frequency range | 10uHz~50kHz |
| Frequency accuracy | <0.01% |
| Frequency resolution | 5000/decade |
| Amplitude | 0.5mV~2.1Vrms (potentiostatic) 0.1~70% f.s.(Galvanostatic) |
| Mode | Static EIS: Potentiostatic, Galvanostatic, Pseudogalvanostatic, OCP Dynamic EIS: Potentiodynamic, Galvanodynamic Fixed frequency impedance: Potentiostatic, Galvanostatic, Potentiodynamic, Galvanodynamic Multisine EIS: Potentiostatic, Galvanostatic Intermittent PEIS/GEIS |

| Interfaces for System | |
|--------------------------|--|
| Auxiliary port | |
| Digital output | 3(open collector) |
| Digital input | 1(photo coupler) |
| Auxiliary voltage inputs | 3 analog inputs: ±10V For measurement of WE vs. CE CE vs. RE or other signal |
| Analog output | 1 analog output: ±10V For stirrer, MFC, RDE, etc. |
| Misc. port | |
| Peripheral communication | I2C to control external devices |
| Temp. measurement | 1 K-type thermocouple input |

| Software | |
|--------------------------|---|
| Max. step per experiment | 1000 |
| Shutdown safety limits | Voltage, current, temperature, etc. |
| Max. sampling rate | 2usec or 3usec depending on data point number |
| Min. sampling time | Unlimited |
| Sampling condition | Time, dv/dt, dl/dt, temperature, etc. |

| PC Requirement | |
|------------------|--------------------------------|
| Operating system | Windows 7/8/10(32bit/64bit OS) |
| PC specification | Pentium4, RAM 1GB or higher |
| Display | 1600x900 high color or higher |
| USB | High speed 2.0 |

| General | |
|------------------------|----------------------------------|
| Dummy cell | One external dummy cell included |
| Thermocouple | K-type, 1.5 meter long(option) |
| Auxiliary cable | Option |
| Misc. cable | Option |
| Impedance analysis S/W | ZMAN™ software |
| DC data analysis S/W | IVMAN™ software package |

The specifications are subject to change without notice.
Windows is a registered trade mark of Microsoft Corporation.


Designed by

ZIVE LAB
www.zivelab.com



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-2-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified