

High Voltage Electrochemical Workstation SHP480HV



For High voltage Energy Storage &
Conversion Devices Application
-1V~+65V @ 5A
10uHz to 500kHz EIS capability

*Battery pack
Fuel Cell stack
Supercapacitors
Solar Cells*

High Voltage Electrochemical Workstation SHP480HV

The **SHP480HV series** are for high voltage energy storage and conversion devices application.

These are power potentiostat/galvanostat/impedance analyzer which is the best choice complete DC and impedance characterization of high voltage devices such as battery pack, fuel cell stack, supercapacitors and solar cells, etc.

The system is designed under FPGA and DSP control with high speed capability.

DAC Control

: Two sets (for offset & scanning) of high speed 16bit DAC(50MHz) and one set of 16 bit DAC(1MHz) for auxiliary analog output control.

ADC Reading

: Two sets of 16 bit 500kHz ADC for reading voltage/current and 4 channel 16 bit 250kHz ADCs for auxiliary data input such as temperature, auxiliary voltage, etc. It provides high frequency EIS, fast pulse techniques and high speed sampling time.

The **SHP480HV series** provides impedance measurements over the frequency range 10uHz to 500kHz. These systems major target application is for unipolar high voltage application under 480Watt.

With various advanced software packages, user can widen **SHP480HV series** flexibility.


System Features

- Versatile all in one power Potentiostat/Galvanostat/Impedance analyzer. (Not requiring external booster)
- 16 EIS techniques capability including multisine & real time EIS at affordable price
- IR measurement (fixed frequency sine wave) during charge/discharge control is available
- Main applications
 - Batteries, Fuel cell, Supercapacitors, Solar cells, etc.
- User defined alias and unit display for max 3 auxiliary signals simultaneously
- High speed data sampling time
 - 2usec or 3usec depending on data point number
- 3 measurement/control voltage ranges
- 7 measurement/control current ranges
- Internal 542,000 data point storage and continuing experiment regardless of PC failure.
- Full software packages are included as standard
 - EIS test software package(EIS)
 - Energy software package(BAT)
 - Corrosion test software package(COR)
 - Electrochemical analysis software package(EAS)
- Free software upgrade

Hardware Features

- 0V to 65V@5Amp
- Compact design in high power hardware.
- Independent operation by DSP with FPGA
- Simultaneous 3 auxiliary voltage measurements
- Temperature measurement as standard
- 1 auxiliary analog output
- 2 digital outputs & 1 digital inputs
- Separated power and sensing line
- Max. 200kHz EIS capability
- Emergency stop button
- CATIII 1000Volt cell cable with standard 1000Volt alligator clip


• Front View



• Rear View



• Smart LCD Display



Versatility

The **SHP480HV series** come with additional 3 analog inputs (auxiliary voltage input) and 1 analog output along with 2 digital outputs and 1 digital input, and one temperature input for K type thermocouple. It will help users expanding usage of the instrument.

For example,

1. User can measure the additional 3 voltage input by using auxiliary voltage input. (User can measure the voltage between reference and counter electrode and between working and counter electrode as well.)
2. With analog output, the system can control external devices (MFC flow rate, etc.) by $\pm 10V$ full scale.
3. User can control using 3bit Digital Output and External control using 1 bit Digital Input.

Safety and Maintenance

1. Even though the communication failure occurs between PC and SHP480HV, the system continues its experiment on channel and saves the data in SHP480HV memory up to 542,000 data point set. After the communication is restored, SHP480HV will transfer saved data set to PC automatically or user can transfer data set when he/she wants. This function will be highly efficient for long-term experiment.
2. User can define a safety condition setting by inputting his/her own safety levels for voltage, current, temperature etc. If the measurement value exceeds this setting value, the system will automatically stop to protect the system and cell.
3. If the control value of voltage or current is different from measured value, the experiment will stop automatically to protect the cell.
4. Automatic calibration function is available for user calibration.
5. The system is controlled from a PC via USB.
6. There is an emergency button to cell off for emergency.

Application

The SHP480HV power electrochemical workstation is the best choice to complete DC and impedance characterization of high voltage energy storage & conversion devices such as high capacity battery pack, fuel cell stack, supercapacitors and solar cells, etc.

Batteries



SHP480HV model is very well adapted for high voltage battery pack application. It provides various control modes for battery cycling with max 5 Ampere @65Volts. It can support EVS (electrochemical voltage spectroscopy)/GITT/PITT test. Fast pulse capability for GSM, CDMA test is included in battery test software package. Pulse profile measurement function to check pulse shape is available. For ripple simulation test, sine wave charging/discharging with profile display is available.

Fuel Cells




SHP480HV is ideal for characterizing the fuel cell stack and anodic/cathodic process mechanism at development and research grade. This system can be directly used for PEMFC, DMFC, and DEFC etc. I-V curve measurements in a full range of available current (autorange option is active during the I-V scan in order to ensure measurement with continuously high resolution).

Supercapacitors



SHP480HV model has fast potentiostat circuit with high speed data acquisition (50usec/point, burst mode). Also it provides high current EIS test for high voltage impedance sample. These functions are well applicable to super capacitor testing. Charging/discharging capability is used for this application.


Solar Cells



Solar cell development and production requires extensive material and device testing to improve efficiency and match individual cells for panel construction. The SHP480HV is the best solution for photovoltaic cell characterization. With system's AI, AO, DI, and DO, the system can monitor other device's signal and also control them.

Impedance(ESR) measurement for QC/QA

The SHP480HV is also suitable for impedance measurement at fixed frequency using independent software.




You can do the above function by replacing file name. It provides simple results with text data files showing "pass"/"fail" decision. The measurement operation is simply click the "start" button or press space bar.

Main Software

The Smart Manager (SM) is to control SHP480HV model and it provides user defined sequential test by using sequence file, technique menu and batch file. The batch file allows the users to do a serial test by combining sequence files and/or technique files. The SM software is easy to use and supports various electrochemical experiments including functions of system control, schedule file editor, real time graph, analysis graph, user calibration, and data file treatment, etc.


High Voltage Electrochemical Workstation SHP480HV



Technique list

■ Sequence editor

User can design his/her experiment procedure by using TASK sequential routine editor.



Sequence editor

- Constant potential, current, C-rate, power, load, OCP
- Sweep potential, current
- Fast sweep potential, current
- Staircase potential, current
- CC-CV, CP-CV, CL-CV, Crate-CV control
- Id, Is control
- EIS control
- Pulse or sinewave control
- Rest(voltage monitoring only)
- Loop(cycle) control

• Cut-off(Vertex) Condition

- Time(step, test, loop, cycle)
- Current, current density
- Voltage
- Capacity
- C-rate
- -dV
- |dV/dt|
- |dI/dt|
- Aux1
- Eoc
- etc.

Item	OP	Delta/value
Step Time	>=	30
None		
Step End		
Step Time		
Current		
I Density		
Voltage		
[Capacity]		
-dV		
dV/dt		
dI/dt		
Temp. (C)		
AUX1		
AUX2		
AUX3		
Test Time		
Loop Time		
Cycle Time		
Eoc		
[Whr]		
LCC(%)		
LCD(%)		
FCC(%)		
FCD(%)		
Power(W)		
SumQ(Ah)		
SumE(Whr)		
Loop Next		
DI Ch		
Crate		


Cut-off condition

• Sampling Condition

- time, |dI/dt|, |dV/dt|, |dT/dt|, |dA1/dt|, burst time

• Flow View

- User can see the sequence flow at a glance.




Flow view


• Control Task Parameters

Control Mode		
constant	GSTAT	constant current control
	Crate	constant Crate control
	PSTAT	constant voltage control
	POWER	constant power control
	LOAD	constant load control
	CC-CV	constant current constant voltage control
	Crate-CV	Crate constant voltage control
	CP-CV	constant power constant voltage control
	CL-CV	constant load constant voltage control
	Id	Id control
Step	Is	Is control
	OCP	OCP control
	GSTAT	current step control
	PSTAT	potential step control
Sweep	GSTAT	current sweep control
	FAST-G	fast current sweep control
	PSTAT	potential sweep control
	FAST-P	fast potential sweep control
	ACV	AC voltammetry
EIS	GSTAT	galvanostatic EIS
	PSTAT	potentiostatic EIS
	OCP	OCP EIS
	PSEUDO	pseudo galvanostatic EIS
	HFR G	galvanostatic HFR
	HFR P	potentiostatic HFR
	MisineG	galvanostatic multisine EIS
	MisineP	potentiostatic multisine EIS
	RTIG	galvanostatic real time EIS
	RTIP	potentiostatic real time EIS
Pulse	Vpulse	voltage pulse control
	Ipulse	current pulse control
	GSINE	current sine wave control
	PSINE	potential sine wave control
	Rest	rest control
ZRA		ZRA control
Loop		loop control
Device	D OUT	Digital output control
	A OUT	Analog output control
	TEMP CTRL	Temperature controller control
	MUX	Multiplexer control


• Battery charging test using sequence editor



Sinusoidal ripple charging



Sine wave profile




Smart Manager Techniques Software Package

For a wide range of application, techniques software packages for specific experimental techniques are available as standard. Each software package's upgrade will be provided at free of charge.

Basic Techniques

Basic techniques with standard functions


- 1) Potentiostatic
- 2) Galvanostatic
- 3) Double step potentiostatic
- 4) Double step galvanostatic
- 5) OCP measurement
- 6) Potential sweep
- 7) Current sweep
- 8) Cyclic voltammetry
- 9) Fast potential sweep
- 10) Potentiostatic Ru measurement
- 11) Galvanostatic Ru measurement



The above functions can be used sequentially by step control function.


EIS Software Package(EIS)

1. Potentiostatic EIS
2. Galvanostatic EIS
3. Pseudo galvanostatic EIS
4. OCP^(*) EIS
5. Potentiodynamic PEIS
6. Galvanodynamic GEIS
7. Potentiodynamic HFR
8. Galvanodynamic HFR
9. Potentiostatic HFR monitor
10. Galvanostatic HFR monitor
11. Multisine potentiostatic EIS
12. Multisine galvanostatic EIS
13. Intermittent potentiostatic EIS
14. Intermittent galvanostatic EIS
15. Real time potentiostatic EIS
16. Real time galvanostatic EIS



Potentiostatic EIS data of metal sample soaked in NaCl solution over time, 3D Nyquist plot by ZMAN


(*) The system measures open circuit potential before each frequency change and apply AC sine wave on this potential.




Batch function

User can design batch file including multiple technique files and/or sequence files. With this batch file, user can experiment several techniques/sequence in series automatically.


Index	Enable	Count	Setting Loop	Loop End	On	Schedule File(s)	File Name
1	☒	1	Next	Next		C:\Zive Data\sm\schedule\level1.EIS	
2	☒	1	Next	Next		C:\Zive Data\sm\schedule\cccv.L.CCV	
3	☒	1	Next	Next		C:\Zive Data\sm\schedule\l1.CCV	
4	☒	1	Next	Next		C:\Zive Data\sm\schedule\l2.7v.IPE	
5	☒	1	Next	Next		C:\Zive Data\sm\schedule\l3.IPE	
6	☒	1	Next	Next		C:\Zive Data\sm\schedule\l4.2V.IPE	
7	☒	1	Index-1	Next		C:\Zive Data\sm\schedule\cccv.L.CCV	
8	☒	1	Next	Next		C:\Zive Data\sm\schedule\lcan.CCV	




Galvanostatic EIS




Potentiodynamic PEIS



Rs, Cp & Idc vs Vdc plot



Potentiostatic EIS



Intermittent Potentiostatic EIS


High Voltage Electrochemical Workstation SHP480HV

• Energy Software Package(BAT)

BAT software supports IR measurement.


1. Battery test techniques

- CC/CV test for cycle life test of lithium battery
- CC/CC test for cycle life test of NiCd or NiMH battery
- Discharging test
- EVS(Electrochemical voltage spectroscopy)
- Variable scan rate CV
- Potentiostatic IV curve
- Galvanostatic IV curve
- Steady-state CV
- GITT(Galvanostatic intermittent titration technique) test
- PITT(Potentiostatic intermittent titration technique) test




CC/CV test

CC/CC test




EVS test raw data


EVS graph format
(dQ/dV vs. V)



Variable scan rate CV



GITT Test




PITT Test

- Pulse mode is available for GSM & CDMA profile. Pulse shape profile can be measured by user's demand.



Pulse shape profile monitor (micro seconds order)



Current sine wave (charge ripple simulation)

2. Control mode

- Charge: CC, CC-CV, pulse, sine wave
- Discharge: CC, CP, CR, pulse, sine wave

3. Cutoff condition

- time, voltage, current, power, temperature, auxV etc.

Various battery charge/discharge test is available including pulse discharge for GSM and CDMA application.



• Corrosion Software Package(COR)

Corrosion technique supports IR compensation.



SHP480HV minimum voltage range is limited by minus 1Volt.

1. Tafel(Tafel experiment)
2. Rp(Polarization resistance)
3. Potentiodynamic
4. Galvanodynamic
5. Cyclic polarization
6. Ecorr vs. time
7. Galvanic corrosion
8. RpEc trend
9. Reactivation potential
10. Critical pitting potential
11. Critical pitting temperature(option)
12. ASTM critical pitting temperature(option)
13. Potentiostatic ECN
14. Galvanostatic ECN
15. ZRA mode ECN


Each software package's upgrade will be provided free of charge.



Tafel experiment



Rp
(Polarization
resistance)



CYPOL(Cyclic polarization resistance)

• Electrochemical Analysis Software Package(EAS)

SHP480HV minimum voltage range is limited by minus 1Volt

1. Step techniques

CA(Chronoamperometry), CC(Chronocoulometry), CP(Chronopotentiometry)

2. Sweep techniques

LSV(Linear sweep voltammetry), SDV(Sampled DC voltammetry), Fast CV, Fast LSV

3. Pulsed techniques

DPV(Differential pulse voltammetry), SWV(Square wave voltammetry), DPA(Diff. pulse amperometry), NPV(Normal pulsed voltammetry), RNPV(Reverse normal pulse voltammetry), DNPV(Differential normal pulse voltammetry)


Control & Real Time Graph

Smart Manager provides virtual control panel for control & data acquisition with real time graph.


User can control and monitor for specific channel in detail and he/she can monitor data in VOI(value of interest) window and channel status in one window. Real time graph's X,Y axis format will be changed per technique automatically. It can be defined by user's demand per techniques.

For experiment using sequence file or batch file, user can designate X,Y parameter on three different real time graph. The real time graph's format can be also selected.

The real time graph and VOI will be changed depending on DC test or impedance test automatically. The virtual control panel always displays the graph for recent test result. For impedance measurement, wave monitor will be displayed on real time graph to check wave's quality. This monitor can be switched to Lissajous(I vs. E) plot.



Strip Chart





Strip chart recorder function provides real time graph function independently. You can monitor 2 Y axis data such as voltage, current, auxV1,2,3, temperature, power, and capacity etc. in real time.

■ Graph



Smart Manager's graph function is to simplify the operation. There are 3 kinds of graph per each experiment. You can change X, Y1, Y2, Y3, Y4 axis parameter as you want. Each graph provides shortcut buttons. When you click these buttons, the format of the graph will be changed accordingly.




In DC and Cycle graph, whenever you click or the parameters which are related to current such as current, capacity, energy, power, load, etc., are changed into calculated specific value or density value, respectively.

: value divided by weight

: value divided by active area




1) DC Graph

- For general data display
- 4 shortcut buttons: I vs. V, E vs. I, I vs. time, V vs. Q
- Graph parameters: time, Eref, I, Eoc, Id, Aux1, Aux2, Aux3, temp, LogI, Load, ChQ, DchQ, ChQs, DchQs, Ch P, Dch P, Ch-Wh, Dch-Wh, Sum Wh, Sum Q, Sum |Q|, |Q|, Rp, dQ/dV



2) EIS Graph

- For EIS data display
- 3 shortcut buttons: Nyquist plot, Bode plot, Cs vs. frequency
- Graph parameters: Frequency, Zre, -Zim, Zmag, Zph, Y, Yimg, Y, |Y|, Yph, LogY, LogY, Rs(R-C), Cs(R-C), Rp(R|C), Cp(R|C), Rs(R-L), Ls(R-L), Q(R-L), time, Vdc, Idc, temp, Aux(1,2,3)




3D Bode Plot by ZMAN
Technique used: Potentiodynamic impedance measurement by using a corrosion cell

3) BAT Graph

- For battery cycle data display
- 3 shortcut buttons: cycle capacity, cycle average, Log(cycle No) vs. depth of discharge plot.
- Graph parameters: cycle number, Ch Q, Dch Q, Sum Q, Coulomb Eff, Ch-Wh, Dch-Wh, Sum Wh, Energy Eff, MinV, MaxV, ChQs, DchQ, ChVavg, DchVavg, Vavg

■ Data Export to ASCII & Excel File




Selectable between 'Convert data on graph only' and 'Convert selected file(s)'

Data Analysis Software


ZIVE data file can be used for analysis by using external IVMAN™ software for DC analysis, IVMAN DA™ software for battery data analysis, IVMAN PA™ software for photovoltaic cell data analysis and ZMAN™ software for EIS data analysis without license.

ZMAN™ EIS Data Analysis Software


- Model simulation and fitting
- 2D- and 3D-Bode- and Nyquist plots
- Automatic equivalent circuit model search function
- Project concept to handle multiple EIS data analysis
- Parameter plot from fitted elements value
- Compatible with data format from Zahner, Gamry, Ametek etc. (License code is needed.)
- Various weighting algorithm
- Model library and user model
- KK plot
- Batch fitting for project data
- Impedance parameter simulation
- Interpolate bad data
- Black-Nichols plot
- 3D graph setting option
- Improved model editor
- Application model library for automatic searching
- Parameter simulation of model
- Genetic algorithm option for initial guessing
- Automatic initial guessing
- Trace movie function on fitting
- Free for ZIVE's data format(*.seo, *.wis) analysis (No license code required.)
- Circle fitting
- Data editing available (insert, delete, edit)
- Add/subtract element parameters
- Add/subtract model parameters
- Impedance, Z in polar, admittance, Y in Polar, modulus, M in polar, dielectric constant, E in polar. data display
- Empty cell capacitance calculation
- Find file function
- Data replacement by formula function
- Cursor data display
- Model finding result automatic sorting by Chi square value
- R, C, L, R, Q preview & graphic
- ZHIT function
- Mott-Schottky analysis
- Donor density vs. Vfb graph
- C vs. voltage graph




Parameter plot




Importing 3rd parties ASCII data file




Project manager with data preview




Cursor data display




Data replacement by formula function




Fitting display




Model editor & model library




Automatic model searching




LEVM fitting




Parameter simulation




Element add/subtraction




Finding data file menu



Circular fitting




2D Nyquist plot




2D Bode plot


High Voltage Electrochemical Workstation SHP480HV




Rp,Cp vs frequency (R/C)




Empty cell capacitance




Donor density vs. Vfb graph and analysis




KK consistency




Mott-Schottky analysis window




C/R-V graph




Bode & Nyquist overlay & 3D plots




Measured data




V vs. Q




Cycle graph



dQ/dV vs. V




dV/dQ vs. Q



Export ASCII file



IVMAN™ Photovoltaic Cell Analysis



- Automatic analysis of parameters
- open circuit voltage, open circuit current, max. power, efficiency, photo induced current, diode quality factor, series resistance, etc.

IVMAN™ DC Data Analysis Software



- IVMAN™ software package consists of
- IVMAN software
- IVMAN utilities
- IVMAN main software
- IVMAN differential analysis software
- IVMAN photovoltaic cell analysis
- IVMAN Tafel analysis
- IVMAN extractor
- IVMAN peak find module




IVMAN DA™ Battery Test Data Analysis Software

- Battery test data analysis
- Electrochemical voltage spectroscopy (dQ/dV vs. V)
- Voltage vs. Capacity analysis (V vs. Q)
- Cycle graph (Q vs. cycle)
- Differential voltage graph (dV/dQ vs. Q)




IVMAN™ Main Software

- Ideal for DC corrosion data analysis and electro-analytical data analysis
- Initial guessing function on Tafel analysis
- Polarization resistance fitting
- 3D graph
- Find peak function
- Interpolation, differentiation, integration etc.
- Reporting function



Polarization resistance fitting



Polarization analysis result

High Voltage Electrochemical Workstation SHP480HV

• Battery Jig

Prismatic cell jig



- 4 pin lever type cylindrical cell jig
- 8 channels, Max.current 50A
- Direct connection method for jig cables
- Compatible with 26650 and 21700 batteries
- For high current application

Through-Plane Conductivity Test Jig



- 2-probe type conductivity measurement
- Max. sample thickness: 40 mm
- 4 mm banana plug connection
- For membrane conductivity measurement

• Battery Clamp Cable



• Faraday Cage



- size : 300 x 300 x 398mm (WxDxH)

Specification

Main System	
PC communication	USB2.0 high speed
Line voltage	100~240VAC, 50/60Hz
Power requirement	390Watt
Size	285 x 550 x 207.2mm(WxDxH)
LED indicator	Run, Comm

System	
Cell cable set	1 meter high voltage power cables(WE/CE) & Reference, working sense, AuxV1 cable
Control DAC	DSP with FPGA 2x16bit DAC(50MHz) for bias & scan 1x16bit DAC(1MHz) for analog output
Data acquisition ADC	2x16bit ADCs(500kHz) for voltage, current 4x16bit ADCs(250kHz) for auxiliary voltage and temperature reading
Calibration	Automatic
Filter selection	4ea(5Hz, 1kHz, 500kHz, 8MHz)
Scan rate	0~10V/sec
Internal data memory	542,000 points
LCD display	DC & EIS mode automatically

Power Amplifier(CE)	
Control voltage	-1V to +65V
Max. current	±5A
Control speed selection	4ea
Bandwidth	1MHz
Slew rate	10V/μsec

Potentiostat Mode (voltage control)	
Voltage control	
Control voltage range	-1V to +65V, 0V to +6.5V, 0V to +650mV
Voltage resolution	16 bit per each range
Voltage accuracy	±0.05% f.s.(gain x1)
Max. scan range	-1V~+65V vs. ref. E
Current measurement	
Current range	10 ranges 5A~5nA 5nA with gain
Current resolution	16 bit 150μA, 15uA, 1.5uA, 150nA, 15nA, 1.5nA 150pA, 15pA, 1.5pA, 150fA
Current accuracy	±0.1% f.s.(gainx1)>5uA f.s.

Galvanostat Mode (current control)	
Current control	
Control current range	max. ±5A ± full scale depending on selected range
Current resolution	16 bit 150μA, 15uA, 1.5uA, 150nA, 15nA, 1.5nA 150pA, 15pA, 1.5pA, 150fA
Current accuracy	±0.1% f.s.(gainx1)>5uA f.s.
Voltage measurement	
Voltage range	65V, 6.5V, 650mV
Voltage resolution	16 bit 1.984mV, 198.4uV, 19.84uV
Voltage accuracy	±0.05% f.s.(gain x1)

Electrometer	
Max. input voltage	65V
Input impedance	>1x10 ¹³ Ω 6pF
Bandwidth	>6.5MHz
CMRR	>128dB

EIS Measurement for System	
Frequency range	10uHz~500kHz
Frequency accuracy	<0.01%
Frequency interval setting (Point / Decade)	Max. 1000/decade (<43mHz) Max. 5000/decade (>43mHz)
Amplitude	3.5mV~20Vrms (Potentiostatic) Bias required 0.1~70%pp f.s.(Galvanostatic)
Mode	Static EIS: Potentiostatic, Galvanostatic, Pseudogalvanostatic, OCP Dynamic EIS: Potentiodynamic, Galvanodynamic Fixed frequency impedance: Potentiostatic, Galvanostatic, Potentiodynamic, Galvanodynamic Multisine EIS: Potentiostatic, Galvanostatic Intermittent PEIS/GEIS Real time PEIS/GEIS

Interfaces for System	
Digital & Auxiliary port	
Digital output	3(open collector)
Digital input	1(photo coupler)
Auxiliary voltage inputs	2 analog inputs: ±10V For measurement of WE vs. CE CE vs. RE or other signal
Analog output	1 analog output: ±10V For flow rate control etc.
Peripheral communication	I2C to control external devices
Temperature input port	
Temp. measurement	1 K-type thermocouple input

Software	
Max. step per experiment	1000
Shutdown safety limits	Voltage, Current, Temperature, etc.
Max. sampling rate	2μsec or 3μsec depending on data point number
Min. sampling time	Unlimited
Sampling condition	Time, dv/dt, dl/dt, temperature, etc.

PC Requirement	
Operating system	Windows 7/8/10(32bit/64bit OS)
PC specification	Pentium4, RAM 1GB or higher
Display	1600x900 high color or higher
USB	High speed 2.0

General	
Dummy cell	One external dummy cell included
Thermocouple	K-type, 1.5 meter long(option)
Auxiliary cable	Option
Misc. cable	Option
Impedance analysis S/W	ZMAN™ software
DC data analysis S/W	IVMAN™ software package

The specifications are subject to change without notice.
Windows is a registered trademark of Microsoft Corporation.

Designed by

ZIVE LAB
www.zivelab.com

WonATech

WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-2-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified