



*EChem
Accessories*

ACCESSORIES CATALOG 2021

For Electrochemical Experiments


Contents

For Electrochemical Experiments

● Corrosion Cell Kit	1
● Flat Cell Kit	2
● Plate Cell Kit	3
● Permeation Cell Kit	4
● Photo Echem Cell Kit	4
● Electrode&Electorde Holder	5
Pt Plate Electrode / Universal Electrode Holder / Flat Specimen Holder	
● Faraday cage	6
● Battery Jig	6
Battery Jig / Coin Cell Jig / Pouch Cell Jig / Coin Cell Holder	
● Membrane Conductivity Cell	8
● Through-Plane Conductivity Jig	9
● Fuel Cell Hardware Fixture	10
● Power Booster	10
● Softwares	12
ZMAN / IVMAN / SIM4U	

Corrosion Cell Kit

The CCK series corrosion cell kit is based on a standard glass reaction flask, 1 liter ~ 100ml. All wetted parts are made of chemically resistant materials such as Teflon, Pyrex and SUS 316. The standard cell configuration consists of a cylindrical metal sample working electrode, a graphite rod, a counter electrode, a gas bubbling tube, and a reference electrode in a Luggin capillary. A graphite rod as counter electrode, a reference electrode and a flat specimen holder could be ordered separately as an option. The water-jacketed type corrosion cell kit and alkaline resistance cell kit made with Teflon are also available.



CCK1, Standard Type
Optional FSH2 & Thermometer

WCCK1, Water-Jacketed Type
Optional FSH2 & Thermometer

Specifications

Vial volume (depending on model))	CCK series : 500 ml & 1 liter WCCK series : 500 ml & 1 liter
Cylindrical sample holder material	
Tube	Pyrex®, 6.35 mm dia.
Compression gasket	Teflon®
Cylindrical metal sample	Steel
Chemical compatibility	
Wetted materials	Pyrex®, Teflon®
Non-wetted materials	Above, plus stainless steel and Viton®
Reference electrode(option)	
Type	SCE or Ag/AgCl reference electrode
Size	9 mm diameter OD, 110 mm long
Counter electrode(option)	
Graphite rod	6 mm diameter, 30 cm long
Pt flag counter electrode	active area : Pt plate type : 1cm ² , 4cm ² , 5cm ²
Flat specimen holder(option)	
Specimen diameter	FSH2 : 15.5 mm ~ 22 mm FSH15 : 18.5 mm~25 mm dia.
Specimen thickness	0.3 ~ 5.8 mm

All specifications are subject to change without notice.

Parts Included

For CCK1 & WCCK1

Cell vial	Pyrex®, 1L
Cylindrical metal sample & tube	Steel / Pyrex®, 6.35mm dai. x 4.35mm dia.
Cell clamp	Stainless steel
Luggin capillary	Pyrex®
Gas bubbler	Pyrex®
Cell Top	Teflon®
Other miscellaneous parts such as stopper / O-ring	MC Nylon® / Viton®

Please contact us for other replacement parts.

Ordering Guide

Standard type	
1 liter volume	CCK1
500 ml volume	CCK05
Water-jacketed type	
1 liter volume	WCCK1
500 ml volume	WCCK05

Components can vary depending on the type of cells.

Optional Items

Flat specimen holder	
Active area : 11.28 mm dia.	FSH2
Active area : 15 mm dia.	FSH15
Counter electrode	
Graphite rod, 150mm long	GR002H
Graphite rod, 300mm long	GR002
Reference electrode	
Saturated calomel reference electrode	WA1001
Ag/AgCl reference electrode	WA1004



FSH2

GR002

WA1001

Flat Cell Kit

The flat cell kit was designed to evaluate plate material such as metal(coupons), semi-conducting plate, etc. A sample plate will be placed one sample holder by fixing knob and maximum 300ml sample volume is acceptable. A water jacketed version is also available. A graphite plate which is placed in one side of the cell is supplied with a cell and can be used as a counter electrode. A Luggin capillary is also included while a reference electrode should be purchased separately. Instead of graphite plate, a platinum wire can be also used as counter electrode by putting through either of the ports on cell body. You can select PTC1 or PTC2 plate test cell kit for small solution volume, which is explained on next page.



FCK2 Standard Type



WFCK2 Water-Jacketed Type

Features

- Ideal for testing of flat specimen
- Easy to use
- Fast and easy disassembly
- Detachable counter electrode
- Two opening areas

Applications

- Polarization test
- Galvanic corros
- Electrochemical noise measurement
- EIS measurement
- Cyclic voltammetry

Specifications

Sample test area	
One side	1 cm ²
The other side	5 cm ²
Sample thickness	Up to 20 mm
Cell volume	up to 300 ml
Material	
Cell body	Pyrex®
Cell end	Polycarbonate
O-ring	Viton®

All specifications are subject to change without notice.

Ordering Guide

Standard type	FCK2
Water-jacketed type	WFCK2

Optional Items

Reference electrode	
Saturated calomel reference electrode	WA1001
Ag/AgCl reference electrode	WA1004

Parts Included

Graphite plate	GR001
Luggin capillary	LGFC

Please contact us for other replacement parts.

Plate Cell Kit

The plate test cell kit, PTC1, is designed to evaluate plate material such as metal(coupons), semi-conducting plate, etc. In evaluation, a sample plate will be placed between two cell blocks. A counter electrode (graphite rod or Pt wire type) and a reference electrode should be ordered separately.

- PTC1 has a electrode holder part, a solution block part, a bottom block part and a thickness adjustment dial knob.
- The active area, which is to be exposed to electrolyte, can be selected by O-ring's position.



PTC1



Active area will be determined by O-ring's position.

Specifications

Sample test area	width: >15mm, thickness: 0.1~10mm
Materials	Teflon®
Active area	
Using small O-ring	1 cm ²
Using large O-ring	5 cm ²

All specifications are subject to change without notice.

The plate test cell, PTC2, is a simple cell for electrochemical testing of coated samples. Also it can be a perfect choice for measuring EIS(Electrochemical Impedance Spectroscopy) of painted metal specimens. The PTC2 is very easy to assemble.



PTC2

Specifications

Sample	
Size	60x60mm or more
Thickness	>7mm
Dimensions	
Base	Approx. 132x90x10mm(WxDxH)
Cell body	
- internal diameter & length	31.5mm, 80mm long
Hole diameter	9.3mm dia. & 6.5mm dia.

All specifications are subject to change without notice.

Parts Included

Cell body	Pyrex®
Base and cell top	Teflon®
Cell clamp	Stainless steel
O-ring	Viton®

Ordering Guide

Plate test cell	PTC1
-----------------	------

Ordering Guide

Plate test cell	PTC2
-----------------	------

Optional Items For PTC1 & PTC2

Reference electrode	
Saturated calomel reference electrode	
- 9mm OD, KT glass tip	WA1001
Ag/AgCl reference electrode	
- 9mm OD, KT glass tip	WA1004
Counter electrode	
Graphite rod	
- 6mm dia. 15cm long	GR002H

Permeation Cell Kit

The permeation cell kit, PMC1, is a spinoff of flat specimen cell kit, FCK2 series, and is designed for permeation test. A membrane or a permeation foil can be placed between two glass half cells.

Two graphite plates which can be used as counter electrode and two Luggin capillary are included as standard. Membrane and reference electrode should be ordered separately.


Specifications

Sample test area	
One side	1 cm ²
The other side	5 cm ²
Dimensions	
Cell vial volume	150ml x 2 ea
Chemical Compatibility	
Wetted materials	Pyrex®, Polycarbonate

All specifications are subject to change without notice.

Ordering Guide

Permeation Cell Kit - Standard type	PMC1
Permeation Cell Kit - Water-jacketed type	WPMC1



Permeation Cell Kit, Standard Type



Permeation Cell Kit, Water-Jacketed Type

Photo Echem Cell Kit

The photoelectrochemical cell having a wide optical window is designed to characterize electrode material under lighting condition. The 2 or 3 electrode test is available. Based on a standard model, PCELL1, the attachments are interchangeable between cells according to user's applications. It is a gas tight sealed cell.



PCELL1
- Standard Model

Specifications


Materials	cell body: PEEK optical window: quartz glass others: SUS 304, Viton O-ring	
Dimensions	optical window dia. 18mm cell dimensions 74.3x40x110mm(WxDxH) (PCELL1)	
Electrolyte volume	max. 6ml (PCELL1)	
Sample size	for PCELL1&2 width: >25mm height: 25~62mm	for PCELL3 width: <18mm height: <22mm
Counter electrode	coiled Pt wire (included)	
Reference electrode	6mm OD electrode available (option)	

All specifications are subject to change without notice.

Ordering Guide


PCELL1 - Standard

- standard type
- one optical window mounted in front of electrolyte chamber



PCELL2

- two optical windows arranged to face each other
- suitable for absorbance measurement with a transparent electrode




PCELL3

- cell kit with a specimen holder
- small sample can be fixed inside the electrolyte chamber



PCELL4

- cell kit with a specimen holder
- small sample can be fixed inside the electrolyte chamber



Electrode & Electrode Holder

Pt Plate Electrode

The Pt plate electrode is made up with a supporting rod with a suitable sized piece of Pt foil at the end. The rod is shielded by a glass tube. This Pt electrode can be served as working electrode or counter electrode.

Specifications

Pt foil	
Active area	1cm ² , 4cm ² , 5cm ²
Pt plate thickness	0.2mm
Rod	
Material	stainless steel
Isolated glass tube	6mm dia.
Length	250mm(rod + contact pin)

All specifications are subject to change without notice.



Ordering Guide

Active area	Pt plate
1cm ²	PFL1
4cm ²	PFL4
5cm ²	PFL5

Universal Electrode Holder

The universal electrode holder, UEH1, is designed to hold various sizes of electrode. The UEH1 has 4 holes to hold electrodes and three of them have a screw to adjust its hole size. The hole size is available from 1.6mm to 10mm. The material of plate is Teflon®, which has high resistance to chemicals and its white color helps user to recognize a tiny change of samples during experiments.

Specifications

Holes	
Number of holes	4
Hole size	1.6mm dia. x 1ea 6.2mm dia. x 1ea 9.6mm dia. x 1ea 10mm dia. x 1ea
Rod	
Material	stainless steel
Diameter	6mm diameter
Length	Max. 150mm

All specifications are subject to change without notice.



Universal electrode holder, UEH1, with optional electrodes and glass vial

Ordering Guide

Universal electrode holder	UEH1

Flat Specimen Holder

The FSH series are sample holders to accommodate flat specimens.

- Pyrex® tube : 6.3mm dia.



Ordering Guide

Flat Specimen Holder Active area : 11.28mm dia. Sample size : 15.5mm~22mm dia. / 0.3~5.8mm thickness	FSH2
Flat Specimen Holder Active area : 15mm dia. Sample size : 18.5mm~25mm dia. / 0.3~5.8mm thickness	FSH15

Faraday Cage

The faraday case, Farad2, is an essential item for electroanalytical experiments. It is well designed to block out external EMI noise and firmly enclosure all the components of electrochemical cell (electrodes, vials, etc.). The spacious interior allows you to set up electronic components or systems easily.

Specifications

Material	
Exterior	powder-coated steel
Interior	powder-coated steel with Teflon®-coated bottom
Window	fine SUS mesh embedded in acryl plates
Access	
Number of holes	2
Size	30mm dia.
Position	right hand side and back side
Dimensions	
Overall	318 x 311 x 409mm(WxDxH)
Window	100x300mm(WxH)

All specifications are subject to change without notice.



Ordering Guide

Faraday cage	Farad2
--------------	--------

Battery Jig

Battery & Coin Cell Jig

- Easy to hold cylindrical cell and/or coin cell
- Wide contact point with gold coated contact area
- 4 contact point type(Kelvin probe) is available to minimize voltage drop or 2 contact probe type.
- Individual channel operation is available.
- Rack type is available.



UCJ1

- Jig for a single coin cell
- Easy to hold a coin cell by pulling the lever
- Banana connector type



CCJ8F2-8PS


- Jig for 8 channel coin cells
- 2 probe type coin cell jig for WBCS series battery cycler
- 8ch per panel
- S type cell connector

Ordering Guide

	Part No.	Description
1st	CCJ	Coin cell jig* ¹⁾
	UCJ	Universal cell jig* ²⁾
2nd	Channel No.	Number of channels required
3rd	F2	2 pin probe
	F4	4 pin probe
4th	Z	4 cell banana connector
	H	For High temperature
	-	Normal type
5th	-L	Lever type
	-	Normal type
6th	4P	4ch per panel
	8P	8ch per panel
7th	S	S type cell connector
	L	L type cell connector
	M	M type cell connector

1) Height fixed type only 2) Height adjustable

Battery Jig



CCJ8FZ2-8P

- Jig for 8 channel coin cells
- 2 probe type coin cell jig for ZIVE series, general battery cycler and potentiostat having banana connectors
- 8ch per panel
- Gap between the pins cannot be adjustable.



CCJ200F2-L20P


- Jig for 200 channel coin cells
- 2 probe type coin cell jig for WBCS series battery cycler
- Lever type
- 20ch per panel

Pouch Cell Jig

- ⊙ Pull-down contact type
- ⊙ 4 contact point type(Kelvin probe)



PPJ4F4-4PS



PPJ1(1ch)



PPJ24-4P

Ordering Guide

Contact Type	Part No.
Pull-down contact type	PPJ*1)

1) * : number of channels

Coin Cell Holder

For WPG/WMPG/WBCS System

- ⊙ Direct connect to cell connector



CCH2L



CCH2

For ZIVE System

- ⊙ D-SUB connector type



Ordering Guide

Description	Part No.
For low current model - WMPG1000Ls/Le/Lx, WBCS3000Ls/Le/Lx series	CCH2L
For standard current model - WPG, WMPG, & WBCS3000S series	CCH2

Ordering Guide



Description	Part No.
For CR2032 coin cell	CCH3-20
For CR2450 coin cell	CCH3-24

Membrane Conductivity Cell


The membrane conductivity cell, MCC, is designed to measure ionic conductivity by simply loading a membrane into cell hardware. The MCC adopts 4 point probe for measuring conductivity. By passing current through two outer electrodes and measuring the voltage through the inner electrodes, it allows the measurement of the conductivity. In the 4-electrode configuration, there is virtually no current flow at the inner voltage sensing electrodes. Therefore, polarization does not occur. The second benefit of the 4-electrode sensor is its tolerance of electrode coating. Since the 4-electrode technique measures potential drop rather than resistance, the measurement remains accurate, despite minor coating. The 2 probe measurement is also available by attaching the working and sensing electrical connections to the cathode side while attaching the counter and reference electrical connections to the anode side. Please see the below configuration.

By placing the conductivity cell between the anode and cathode conduction plate, you can simply assemble the conductivity cell into your fuel cell hardware.

- Supports 2 or 4 electrode measurement
- Material
 - cell body : PEEK
 - wire : platinum
- Operating temperature : to 130 °C
- Fuel cell hardware available
 - : 5, 25 cm² fuel cell test hardware (not included, provided by WonATech)
- Easy to assemble



Connecting for a 2-electrode measurement



Connecting for a 4-electrode measurement

Specifications

Material	
Cell, clamp & nut	PEEK
Electrode(S-wire/Z-wire)	Platinum
Dimensions	
Conductivity cell	76.2x76.2x20 mm(WxHxD)
Conductivity clamp	48x50x7 mm(WxHxD)
S-wire (inner electrodes)	84 mm long x 1.0 mm dia.
Z-wire (outer electrodes)	120 mm long x 1.0 mm dia.
Access	
Voltage measurement (S-wire)	two, inner ports
Current measurement (Z-wire)	two, outer ports
Temperature measurement	one, side port

All specifications are subject to change without notice.

Ordering Guide

Membrane conductivity cell	MCC
----------------------------	-----

Through-Plane Conductivity Jig

Because the conductivity of a material is directly linked with ohmic losses, the measurement of ionic conductivity is crucially important in order to evaluate the performance of a newly synthesized material such as ion exchange membrane(IEM) and proton exchange membrane(polymer electrolyte membrane, PEM).

Today ion exchange membranes are receiving considerable attention and are successfully applied for desalination of sea and brackish water and for treating industrial effluents. And proton exchange membrane(PEM) is one of the key components for various consumer related applications for fuel cells, e.g. automobiles, back-up power, portable power etc. For example, in PEMs, protons can transport in two directions, across the membrane and through the membrane. This results in two conductivities, in-plane conductivity and through-plane conductivity. For PEM fuel cells, through-plane conductivity measurement is more meaningful than in-plane because proton transfer occurs in the through-plane direction.


The conductivity of the membrane can be calculated based on the measured resistance by the following equation:

$$\sigma = \frac{L}{RWT}$$

where σ is the membrane conductivity(S/cm), L is the length between the electrodes, R is the measured resistance, W is the membrane width, and T is the membrane thickness.



MCJ1 (Through-Plane Conductivity Test Jig)



The MCJ1 Through-plane conductivity test jig helps user to setup a 2-probe electrochemical cell consisting of 2 stainless steel probes that sandwiches the membrane to measure through-plane conductivity of membranes. The MCJ1 is designed to hold a membrane by pulling a lever.

Normally, a number of galvanostatic alternating current(AC) electrochemical impedance spectroscopic (EIS) techniques or DC techniques are used to estimate the membrane conductivity. User can set up a perfect system with one of ZIVE series Electrochemical Workstation with MCJ1 conductivity test jig for through-plane conductivity measurements.

Specifications

Sample size	>30mm diameter
Sample thickness	max. 40mm
Sample contact material	304 stainless steel
Overall dimensions	70 x 135 x 174mm(WxDxH)
Connection	4mm banana plug

All specifications are subject to change without notice.

Ordering Guide


Through-plane conductivity jig	MCJ1
--------------------------------	------

Fuel Cell Hardware Fixture

- For PEMFC, DMFC
- Max.operating temperature : 120 °C or 180 °C
- Active area : 5, 25cm²
- Components : serpentine flow pattern, cartridge heater(2), current collector(2), cell graphite(2), end plate(2), connector
- Thermal Jacket is available as an option.
- MEA is not included.

Ordering Guide


Active Area	Max. Temp.	Part No.
5cm ²	120 °C	SCFC5
25cm ²	120 °C	SCFC25
5cm ²	180 °C	SCFC5H
25cm ²	180 °C	SCFC25H




• SCFC*** fuel cell hardware



• SCFC***H fuel cell hardware



• single serpentine for SCFC5



• triple serpentine for SCFC25

Power Booster

The ZIVE ZB series boosters are a new generation of single or multi-channel high current instrumentation and they are designed to increase the maximum current and/or maximum voltage of ZIVE series potentiostat/galvanostat.

The ZIVE ZB series boosters have full dc capabilities and are ideal for a wide range of electrochemical applications including high speed voltage/current pulse techniques. And impedance analysis techniques such as single- and multi-sine and HFR test, etc. are also available. Wide frequency ranges covering 10uHz to 1kHz(10kHz) depending on system power enables user to characterize energy storage devices and electrochemical cells over their full frequency range.

This ZIVE ZB series boosters are designed as stand alone type or rack mounted type and have multiple booster modules placed inside it. The power capability can be growing by adding module units to the existing system (factory configuration).



ZB1

ZB2

ZB3

- For high voltage/high current application
- Modular type design
- EIS capability
- Sine wave simulation available
- Simple operation and accurate result
- Safety features for user and instrument itself



ZB4



Rack Type

Power Booster

Specification

Control & Measurement	
Maximum Power Dissipation	3,840Watt
Maximum Current	200A
Minimum Frequency	10uHz
Maximum Frequency	1kHz ~ 10kHz (depending on power)
Current Range	single
Voltage Range	single
Input Impedance	10 ¹³ Ohm
Accuracy	0.05% ~ 0.1% f.s. (depending on power)
Resolution	16 bit
Rise Time	5usec ~ 500usec (depending on power)
Cooling Method	forced air flow
Data Acquisition	>50usec

* T* This booster needs ZIVE workstation, booster interface cable and cell cable.
The specifications are subject to change without notice.

Ordering Guide


Housing	Model	Max.V	Max.I(<-2V)	Max. I (Bipolar)	Power Dissipation(Watt)
ZB1	ZB530B	5V		30A	450
	ZB1030U/1020B	10V	30A	20A	459/480
	ZB2015U/2010B	20V	15A	10A	409/480
	ZB408U/405B	40V	9A	5A	410/480
ZB2	ZB560B	5V		60A	900
	ZB1060U/1040B	10V	60A	40A	918/960
	ZB2035U/2020B	20V	35A	20A	955/960
	ZB4015U/4010B	40V	15A	10A	770/960
ZB3	ZB1090U/1060B	10V	90A	60A	1,377/1,440
	ZB2050U/2030B	20V	50A	30A	1,365/1,440
	ZB4025U/4015B	5V	25A	15A	1,365/1,440
ZB4	ZB1080B	10V		80A	1,920
	ZB2060U/2040B	20V	60A	40A	1,683/1,920
	ZB4030U/4020B	40V	30A	20A	1,539/1,920
ZBR2	ZB5190B	5V		190A	3,800
	ZB10160B	10V		160A	3,840
	ZB20120U/2080B	20V	120A	80A	3,480/3,840
	ZB3090U/3030B	30V	90A	30A	3,447/2,160
	ZB4070U/4035B	40V	70A	35A	3,591/3,360

Model Name ****B is for voltage bipolar type, ****U is for voltage unipolar type [minimum voltage -1V or -2V(ZB20120U)]

* Customized specification is available. Please contact WonATech sales team.

EIS Data Analysis Software, ZMAN™

- Model simulation and fitting
- 2D- and 3D-Bode- and Nyquist plots
- Automatic equivalent circuit model search function
- Project concept to handle multiple EIS data analysis
- parameter plot from fitted elements value
- compatible with data format from Zahner, Gamry, Ametek etc. (License code is needed)
- Various weighting algorithm
- Model library and user model
- KK plot
- Batch fitting for project data
- Impedance parameter simulation
- Interpolate bad data
- Black-Nichols plot
- 3D graph setting option
- Improved model editor
- Application model library for automatic searching
- Parameter simulation of model
- Genetic algorithm option for initial guessing
- Automatic initial guessing
- Trace movie function on fitting
- Free for ZIVE's data format(*.seo, *.wis) analysis (no license code required)
- Circle fitting
- Data editing available (insert, delete, edit)
- Add/subtract element parameters
- Add/subtract model parameters
- Impedance, Z in polar, admittance, Y in Polar, modulus, M in polar, dielectric constant, E in polar. data display
- Empty cell capacitance calculation
- Find file function
- Data replacement by formula function
- Cursor data display
- Model finding result automatic sorting by Chi square value
- R, C R, L R, Q preview & graphic
- ZHIT function
- Mott-Schottky analysis
- Donor density vs. Vfb graph
- C vs. voltage graph



Software

DC Data Analysis Software, IVMAN™


IVMAN™ software package consists of

- IVMAN software
- IVMAN utilities
 - IVMAN main software
 - IVMAN differential analysis software
 - IVMAN photo voltaic cell analysis.
 - IVMAN Tafel analysis
 - IVMAN extractor
 - IVMAN peak find module



IVMAN TA™ Tafel Analysis

- Simple Tafel calculation




IVMAN™ Main Software

- Ideal for DC corrosion data analysis and electro-analytical data analysis
- Initial guessing function on Tafel analysis etc.




IVMAN DA™ Battery Test Data Analysis Software

- Battery test data analysis
- Electrochemical voltage spectroscopy (dQ/dV vs. V)
- Voltage vs. Capacity analysis (V vs. Q)
- Cycle graph (Q vs. cycle)
- Differential voltage graph(dV/dQ vs. Q)




Measured data

V vs. Q




Cycle graph

dQ/dV vs. V




dV/dQ vs. Q

Export ASCII file




Polarization resistance fitting

Polarization analysis result



Time graph

Find peak menu



CV graph


3D graph

IVMAN EX™ Extractor

- Extracting data by cycle number or step
- Exporting ASCII file




IVMAN™ Photovoltaic Cell Analysis




IVMAN PF™ Peak Find Module

- Independent peak finding software



Simulation Software for Cyclic Voltammetry, SIM4U Freeware

- Single or multiple charge transfer steps and first and second-order chemical steps can be used
- Cyclic voltammetry method is used for simulation
- 1D simulation of semi-infinite diffusion processes is used
- The pre-equilibrium can be applied before simulation
- The effect of uncompensated resistance and double layer capacitance can be simulated.
- Measured data and simulated data can be seen together in the plot



Designing the Solution for Electrochemistry



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com
www.zivelab.com

Local Distributor



ISO 9000 & ISO 14000 Qualified