




2 pin pull down type Coin Cell Jig

Pull down type 2 Pins Coin Cell Jig

- Easy to hold coin cell
- Wide contact point with noble coated contact area
- Individual channel operation is available.
- Max. height of coin cell: 6mm
- Max. current: 2A

ex)



Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PP	2 Pin PCB type
	4PL	4 Pin Lever type
	4PK	4 Pin Knob type
2nd	CCJ	Coin cell jig
	PCJ	Pouch cell jig
	CBJ	High current cylindrical cell jig
	UCJ	Universal cell jig
	UCBJ	High current universal cylindrical cell jig
	UCJH	High current universal cell jig
	PRCJ	Prismatic cell jig
3rd	Channel No.	Channel Quantity
4th	Z	Banana connector panel
	H	For high temperature
5th		Normal type
	4P	4ch per panel
	8P	8ch per panel
	16P	16ch per panel
6th	20P	20ch per panel
	S	S type cell connector
	L	L type cell connector
	M	M type cell connector



WonATech Co., Ltd.
 7, Neunganmal 1-gil, Seocho-gu,
 Seoul, 06801, Korea
 Phone: +82-2-578-6516
 Fax: +82-576-2635
 e-mail) sales@wonatech.com
 website: www.wonatech.com

Local Distributor

