





# Lever type 4 pins Coin Cell Jig

## Lever type 4 Pins Coin Cell Jig

- Easy to hold coin cell using clamp lever
- Kelvin probe type 4pin contact
- Jig for WBCS3000L/Le/Lx /S and M1/M2 model
- Wide contact area with spring loaded cell
- Banana socket panel is available for 3rd parties cycler or Potentiostat/Galanostat (Option: B)
- For EIS measurement (Option: Z)
- Individual channel operation is available.
- Optional jig cable length: 2m, 3m, 5m, 6m
- Max height of coin cell: 8.5mm
- Max current: 15A (For M type cell connector)  
10A (with optional banana connector)  
5A (EIS measurement jig)



### ● Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PP	2 Pin PCB type
	<b>4PL</b>	<b>4 Pin Lever type</b>
	4PK	4 Pin Knob type
2nd	<b>CCJ</b>	<b>Coin cell jig</b>
	PCJ	Pouch cell jig
	CBJ	High current cylindrical cell jig
	UCJ	Universal cell jig
	UCBJ	High current universal cylindrical cell jig
	UCJH	High current universal cell jig
	PRCJ	Prismatic cell jig
3rd	<b>Channel No.</b>	<b>Channel Quantity</b>
4th	<b>Z</b>	<b>EIS measurement</b>
	B	Banana connector panel
	H	For high temperature
5th	4P	4ch per panel
	<b>8P</b>	<b>8ch per panel</b>
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	<b>L</b>	<b>L type cell connector</b>
	M	M type cell connector



WonATech Co., Ltd.  
7, Neunganmal 1-gil, Seocho-gu,  
Seoul, 06801, Korea  
Phone: +82-2-578-6516  
Fax: +82-576-2635  
e-mail) sales@wonatech.com  
website: www.wonatech.com

Local Distributor

