


Lever type 4 pins Universal Jig

Lever type 4 Pins Universal Jig

- Easy to hold cylindrical or coin cell using clamp lever
- Kelvin probe type 4pin contact
- Jig for WBCS3000S and M1/M2 model
- Wide contact area with spring loaded cell
- Banana socket panel is available for 3rd parties cycler or Potentiostat/Galanostat
- Individual channel operation is available.
- Optional jig cable length: 2m, 3m, 5m, 6m
- Max current: 15A (For M type cell connector)
10A (with optional banana conector)
5A (with optional EIS measurement)
- Height adjustment is available. Max height: 118.8mm

ex)



● Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PP	2 Pin PCB type
	4PL	4 Pin Lever type
	4PK	4 Pin Knob type
2nd	CCJ	Coin cell jig
	PCJ	Pouch cell jig
	CBJ	High current cylindrical cell jig
	UCJ	Universal cell jig
	UCBJ	High current universal cylindrical cell jig
	UCJH	High current universal cell jig
	PRCJ	Prismatic cell jig
	3rd	Channel No.
Z		EIS measurement
4th	B	Banana conector panel
	H	For high temperature
		Normal type
5th	4P	4ch per panel
	8P	8ch per panel
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	L	L type cell connector
	M	M type cell connector



Optional banana connector panel



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

