




Battery Clamp Cable

Battery Clamp Cable

- Can easily fix coin cell, pouch cell, or battery
- Clamp type that connects to terminals in a clip
- Minimizes voltage drop with a Kelvin probe (except CCC clamp)
- Individual channel operation is available

ex)



- Small alligator clip, FA



- Max. current : 5A
- For clamp stand

- Mid alligator clip, UC



- Max. current : 10A
- For clamp stand
- For clamp cable

- Ordering Guide

	Part No.	Description
1st	UCC	Universal clamp cable
	HPCC	Pouch cell clamp cable for high current
	CCC	Coin cell clamp cable
2nd	L/S/M	L/S/M type connector
	MX	Molex type connector (per 4 channels)
	O	O-ring/banana connector combination
	B	Banana connector
3rd	30	Cell cable length(30cm)
	100	Cell cable length(100cm)
	300	Cell cable length(300cm)

- Mid alligator clip, UCC



UCCS100



- UCCS



UCCB30



- UCCM



- UCCB



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
E-mail: sales@wonatech.com
Website: www.wonatech.com

Local Distributor

