


Battery Clamp Stand

Battery Clamp Stand

- Can easily fix coin cell, pouch cell, or battery
- Clamp type that connects to terminals in a clip
- Minimizes voltage drop with a Kelvin probe (except CC clamp)
- Individual channel operation is available

ex)



● Ordering Guide


	Part No.	Description
1st	4PA	4 contact clip type
	2PA	2 contact clip type
	2PA2	Two 2-contact clip type(4 contact)
2nd	FA ^{*1}	Flat type clamp cable
	UC	Universal clamp cable
	CC ^{*1}	Coin cell clamp cable
	PC ^{*1}	Pouch clamp cable
	HPC	High current pouch cell clamp cable
3rd		Channel Quantity
4th	A ^{*2}	Rear panel cable
	B ^{*2}	Front panel cable
5th	8P	8ch per panel
	16P	16ch per panel
	20P	20ch per panel
6th	L/S/M	L/S/M type connector
	MX	Molex type connector (per 4 channels)
	O	O-ring/banana connector combination
	B	Banana connector
7th	30	Cell cable length(30cm)
	100	Cell cable length(100cm)
	300	Cell cable length(300cm)

*1: Cable type supports only banana connectors.


*2 : A and B are only supported for flat type clamp cable(FA).



2PA2- FA64A-8P-MX



Rear panel cable



2PA2- FA32B-8P-MX



front panel cable



WonATech Co., Ltd.
7 Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Tel: +82-2-578-6516
Fax: +82-2-576-2635
E-mail: sales@wonatech.com
Website: www.wonatech.com

Local Distributor



ISO 9000 & ISO 14000 Qualified