




PCJ1 Single Channel Pouch Cell Jig

Stand alone 4 Pins Pouch Cell Jig

- *Easy to hold pouch cell using clamp lever*
- *Kelvin probe type 4pin contact*
- *Wide contact area with spring loaded cell*
- *Banana connectors panel is for cycler or Potentiostat/Galanostat*
- *Max current: 10A*
Max current: 5A for EIS measurement
- *Terminal width adjustment is available.*



banana socket panel & clamp lever



WonATech Co., Ltd.
7, Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Phone: +82-2-578-6516
Fax: +82-576-2635
e-mail) sales@wonatech.com
website: www.wonatech.com

Local Distributor

