


Knob type 4 pins Prismatic battery/Pouch cell Jig



Knob type 4 Pins Prismatic Battery Pouch cell Jig

- Easy to hold high capacity prismatic battery or pouch cell using knob
- Kelvin probe type 4 pin contact
- Wide contact area with spring loaded cell
- Adjustable cell contact width: 32mm~ 20cm
- High current ring connectors with banana sensing connectors
- Max current: 50A
- Height adjustment by knob is available. Max height: 25cm
- Voltage sense connector: 4mm banana
- Current connector: cable core or ring connector



● Ordering Guide

	Part No.	Description
1st	2PD	2 Pin Pull down type
	2PP	2 Pin PCB type
	4PL	4 Pin Lever type
	4PK	4 Pin Knob type
2nd	CCJ	Coin cell jig
	PCJ	Pouch cell jig
	CBJ	Cylindrical cell jig
	UCJ	Universal cell jig
	UCBJ	High current Cylindrical cell jig
	UCJH	High current Universal cell jig
	PRCJ	Prismatic cell jig
	3rd	Channel No.
4th	Z	Banana connector panel
	H	For High temperature Normal type
5th	4P	4ch per panel
	8P	8ch per panel
	16P	16ch per panel
	20P	20ch per panel
6th	S	S type cell connector
	L	L type cell connector
	M	M type cell connector



WonATech Co., Ltd.
7 Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Tel: +82-2-578-6516
Fax: +82-2-576-2635
E-mail: sales@wonatech.com
Website: www.wonatech.com

Local Distributor

