



UCJ1A Single Channel Universal Jig

Stand alone 4 Pins Universal Jig

- Easy to hold cylindrical or coin cell using clamp lever
- Kelvin probe type 4pin contact
- Wide contact area with spring loaded cell
- Banana connectors panel is for cycler or Potentiostat/Galanostat
- Max current: 10A
Max current: 5A (EIS measurement Jig)
- Height adjustment is available. Max height: 80mm



banana socket panel & height ruler



WonATech Co., Ltd.
7 Neunganmal 1-gil, Seocho-gu,
Seoul, 06801, Korea
Tel: +82-2-578-6516
Fax: +82-2-576-2635
E-mail: sales@wonatech.com
Website: www.wonatech.com

Local Distributor

